





Monoloc 1.5/2.0 Locking Plates System

Surgical Technique

CHANGZHOU KANGHUI MEDICAL INNOVATION CO., LTD

Add. No.11 North Changjiang Road, Xinbei Zone, Changzhou, Jiangsu, China 213022 Tel: 86-519-85195556 Fax: 86-519-85195551 www.kanghui-med.com







Introduction
Indications
Surgical Technique
Preparation
Fracture Reduction and Plate Cont
Screw Insertion
Cortex Screws
Locking Screws
Implant Removal
Tips
Implant
Instruments

Content

Warning:

This description is not sufficient for immediate application of the instrumentation. Instruction by a surgeon experienced in handling this instrumentation is highly recommended.

	P01
	P02

	P03
touring	P04
	P07
	P08
	P11
	P14
	P17
	P18
	P28





Introduction

The Monoloc 1.5/2.0mm Locking Plate system is indicated for fixation of small bones and small bone fragments, including osteotomies, arthrodeses, replantations, reconstructions of small bones and small bone fragments, particularly in osteopenic bone.

All plates have a low profile to reduce softtissue irritation, locking hole design for use with cortex or locking screws, and a cut-to-length feature to minimize inventory.

All locking and cortex screws have a selftapping tip to facilitate insertion and selfretaining StarDrive recess improving torque transmission and increasing resistance to stripping.

Tapered, threaded hole for locking screws

The self-tapping locking screws can be locked in the tapered threaded hole to ensure angular stability. Plate and screw systems in which the screws are locked in the plate function according to the principle of an internal fixator and are used to resolve the following problem situations:

- Primary intraoperative loss of reduction
- Secondary postoperative loss of reduction, particularly in cases of osteoporosis or poor bone quality or of comminuted fractures without bony support
- Compression of the periosteum and the resulting impairment of cortical circulation

Fuctional principle of the internal fixator

When LCP plates are used with angularly stable locking screws, the plate and screws together form a stable system; the stability of the fracture is mainly dependent on the strength of the resulting assembly. Since the plate does not need to be compressed against the bone, blood flow to the bone is not additionally impaired.









Indications

Preparation

1.5 Locking Plates

- Fractures, osteotomies and arthrodeses of
- the phalanges and metacarpals
- Replantations and reconstructions of the phalanges and metacarpals

2.0 Locking Plates

- Fractures of the phalanges, metacarpals,

wrist bones

 Osteotomies and arthrodeses of the interpha langeal joints

This surgical technique describes the application of the LCP Hand system using a Locking Condylar Plate 1.5.

1

Select implant

Select the plates according to the fracture pattern and anatomy of the bone.







Fracture Reduction and Plate Contouring

1

Reduce fracture



Reduce the fracture under image intensification and, if necessary, fix with Kirschner wires or reduction forceps. The reduction method is fracture-specific.

2

Trim plate

Instrument

161210	Cutter
	••••••

Trim the plate to the desired length using bending/cutting pliers. Remove the burrs.

3

Contour plate

Preparation

Instruments

For 1.5 Lo	cking Plates
98092	Bending Pin for LCP Plates 1.5, with thread
For 2.0 Lo	ocking Plates
98081	Bending Pin for LCP Plates 2.0, with thread

If necessary, contour the plate to suit the anatomical conditions. Use bending pins for Locking Plates and thread the pins into the screw holes to contour the plate.

Note:

Do not deform the threaded part of the holes or overbend the plate as this may adversely affect insertion of locking screws. Avoid repetitive bending of the plate.







Preparation

4

Position plate

Instrument

01123 Reduction Forceps with Points

Position the plate over the reduced fracture and, if necessary fix provisionally with Kirschner wires or reduction forceps.



Determine screw type

Depending on the individual case, cortex screws and/or locking screws may be inserted. Determine where locking screws will be used. The locking holes accept both, cortex and locking screws.

If a locking screw is inserted first, ensure that the plate is held securely to the bone to prevent the plate from spinning as the screw locks into the plate.

The final screw placement and the use of locking and cortex screws is determined by the fracture pattern.





Screw Insertion Cortex Screws

1

Pre-drill screw hole for cortex screw

Instruments

For 1.5 Loc	cking Plates
98083	Double Drill Guide 1.5/1.1
98091	Drill Bit Φ 1.1mm, with Quick
	Coupling
98085	Drill Bit Φ 1.5mm, with Quick
	Coupling
For 2.0 Loc	cking Plates
98095	Double Drill Guide 2.0/1.5 for Drill
	Bits 2.0 and 1.5mm
00005	Drill Dit of Emm. with Quick

98085	Drill Bit Ø1.5mm, with Quick
	Coupling
98082	Drill Bit Ф2.0mm, with Quick
	Coupling

Drill the screw hole using the drill bit and the corresponding drill guide / drill sleeve.



2

Determine screw length of cortex screw

Instruments

For 1.5 Locking	Plates
98086	Depth Gauge for Screws 1.5mm,
	measuring range up to 24mm

For 2.0 Locking Plates

98087	Depth Gauge for Screw 2.0 and
	2.4mm, measuring range up to
	40mm

Determine the screw length with the corresponding depth gauge.







Screw Insertion Locking Screws

3

Insert cortex screw

Instruments

For 1.5 Lockin	g Plates
98100	Screwdriver Shaft Stardrive1.5,T4, self-holding, with Quick Coupling
98098	Holding Sleeve 1.5
For 2.0 Lockin	g Plates
For 2.0 Lockin 98099	g Plates Screwdriver Shaft Stardrive 2.0, with Quick Coupling

Insert the self-tapping cortex screw using the corresponding Screwdriver.



1

Pre-drill screw hole for locking screw

Instruments

For 1.5 Lockin	g Plates
98094	LCP Drill Sleeve 1.5, for Drill Bits Φ 1.1 mm
98091	Drill Bit Φ 1.1mm, with Quick Coupling
For 2.0 Lockin	g Plates
98096	LCP Drill Sleeve 2.0, with Scale, for Drill Bits Φ 1.5 mm
98085	Drill Bit Φ 1.5 mm, with Quick Coupling

Drill screw hole through the LCP drill sleeve using the appropriately sized drill bit.insertion of locking screws. Avoid repetitive bending of the plate.







2

Determine screw length of locking screw

Instruments

For 1.5 Locking	g Plates
98086	Depth Gauge for Screws Φ 1.3 to 1.5 mm, measuring range up to 24 mm
For 2.0 Locking	g Plates
98096	LCP Drill Sleeve 2.0, with Scale, for Drill Bits Φ 1.5 mm
98087	Depth Gauge for Screw 2.0 and 2.4 mm, measuring range up to 40 mm

Determine the screw length either with the corresponding depth gauge or with the LCP drill sleeve

When using the LCP drill sleeve check the length directly on the scale of the drill sleeve. Then remove the drill sleeve.



3

Insert locking screw

Instruments

For 1.5 Lockin	ng Plates
98100	Screwdriver Shaft Stardrive 1.5,
	T4, self-holding, with Quick
	Coupling
98098	Holding Sleeve 1.5
899121	Handle, with Mini Quick Coupling
For 2.0 Lockin	ng Plates
For 2.0 Lockin 98099	ng Plates Screwdriver Shaft Stardrive 2.0,
For 2.0 Lockii 98099	ng Plates Screwdriver Shaft Stardrive 2.0, with Quick Coupling
For 2.0 Lockin 98099 98097	ng Plates Screwdriver Shaft Stardrive 2.0, with Quick Coupling Holding Sleeve 2.0
For 2.0 Lockin 98099 98097 899121	ng Plates Screwdriver Shaft Stardrive 2.0, with Quick Coupling Holding Sleeve 2.0 Handle, with Mini Quick Coupling

Insert the self-tapping locking screw using the corresponding screwdriver shaft and the appropriate handle.







Implant Removal

Tips

To remove the plate, first unlock all screws then definitely remove them in a second step. If the screws are not unlocked before removal the plate may rotate while the last screw is being removed and cause soft tissue damage.

Pick-up screws with Stardrive recess without holding sleeve

Select the screw and position the corresponding screwdriver perpendicular to the screw. Pick-up the screw by applying slight axial pressure.

Make sure that the screwdriver remains perpendicular during screw extraction.

Pick-up screws with holding sleeve

All screws can also be picked-up with help of the corresponding holding sleeves.

Slide the holding sleeve onto the screwdriver, until it clicks into place. With the holding sleeve jaws open, mount the appropriate screw onto the screwdriver, then push the holding sleeve until it secures the screw.

Note:

The holding sleeve covers the head of the screw.





Implants

Locking Plate 1.5, straight

Product No.	Size	Material	Model No.	0
29476504	4H	PT	SDXZ40	-
29476506	6H	PT	SDXZ40	

Locking Reconstruction Plate 1.5

Product No.	Size	Material	Model No.	
29476006	6H	PT	SDXZ75	 \sim
29476012	12H	PT	SDXZ75	

Locking Y-Plate 1.5

Product No.	Size	Material	Model No.
29478308	3+8H	PT	SDXZ06

Locking Condylar Plate 1.5

Product No.	Size	Material	Model No.
29475906	2+6H	PT	SDXZ45

Locking Strut Plate 1.5

Product No.	Material	Model No.
24241000	PT	SDXZ68

Locking T-Plate 1.5

Product No.	Size	Material	Model No.	
29476608	3+8H	PT	SDXZ73	
29476408	4+8H	PT	SDXZ74	



333

0

0











Locking Plate 2.0, straight

Product No.	Size	Material	Model No.	
29475404	4H	PT	SDXZ40	
29475405	5H	PT	SDXZ40	
29475406	6H	PT	SDXZ40	
29475407	7H	PT	SDXZ40	
29475408	8H	PT	SDXZ40	

Locking T-Plate 2.0

Product No.	Size	Material	Model No.
29476107	3+7H	PT	SDXZ73



CO O O O O O O O

Locking Condylar Plate 2.0

Product No.	Size	Material	Model No.
29475307	2+7H	PT	SDXZ45



Locking Y-adaption Plate 2.0

Product No.	Size	Material	Model No.
29476307	3+7H	PT	SDXZ06



Locking T-adaption Plate 2.0

Product No.	Size	Material	Model No.
29476207	2+7H	PT	SDXZ72

Rotation Correction Locking Plate 2.0

Product No.	Material	Model No.
24242000	PT	SDXZ78

Mini Condylar Locking Plate 2.0

Product No.	Direction	Material	Model
24243200	Left	PT	SDXZ
24243100	Right	PT	SDXZ

Mini H Locking Plate 2.0

Product No.	Material	Model No.
24244000	PT	SDXZ70











1.5mm Locking Screw Stardrive, self-tapping

Product No.	Size	Material	Model No.
39475806	1.5x6mm	ТА	SDLD01
39475808	1.5x8mm	ТА	SDLD01
39475810	1.5x10mm	ТА	SDLD01
39475812	1.5x12mm	ТА	SDLD01
39475814	1.5x14mm	TA	SDLD01
39475816	1.5x16mm	TA	SDLD01
39475818	1.5x18mm	ТА	SDLD01
39475820	1.5x20mm	TA	SDLD01
39475822	1.5x22mm	TA	SDLD01
39475824	1.5x24mm	TA	SDLD01

1.5mm Cortex Screw Stardrive, self-tapping

Product No.	Size	Material	Model No.	
39475706	1.5x6mm	TA	SSQ01	
39475708	1.5x8mm	TA	SSQ01	
39475710	1.5x10mm	TA	SSQ01	
39475712	1.5x12mm	TA	SSQ01	
39475714	1.5x14mm	TA	SSQ01	
39475716	1.5x16mm	TA	SSQ01	
39475718	1.5x18mm	TA	SSQ01	
39475720	1.5x20mm	TA	SSQ01	
39475722	1.5x22mm	TA	SSQ01	
39475724	1.5x24mm	ТА	SSQ01	

-

Contraction of the second second

2.0mm Locking Screw Stardrive, self-tapping

2.0mm Loci	ang Screw St	ardrive,self-	tapping	
Product No.	Size	Material	Model No.	
39475606	2x6mm	TA	SDLD01	
39475608	2x8mm	TA	SDLD01	
39475610	2x10mm	ТА	SDLD01	
39475612	2x12mm	ТА	SDLD01	
39475614	2x14mm	ТА	SDLD01	
39475616	2x16mm	ТА	SDLD01	
39475618	2x18mm	ТА	SDLD01	
39475620	2x20mm	ТА	SDLD01	
39475622	2x22mm	ТА	SDLD01	
39475624	2x24mm	ТА	SDLD01	
39475626	2x26mm	ТА	SDLD01	
39475628	2x28mm	ТА	SDLD01	-
39475630	2x30mm	TA	SDLD01	-

2.0mm Cortex Screw Stardrive, self-tapping

Product No.	Size	Material	Model No.	
39475506	2x6mm	TA	SSQ01	
39475508	2x8mm	ТА	SSQ01	
39475510	2x10mm	ТА	SSQ01	
39475512	2x12mm	ТА	SSQ01	
39475514	2x14mm	ТА	SSQ01	
39475516	2x16mm	ТА	SSQ01	
39475518	2x18mm	TA	SSQ01	





MonoLoc Locking Plates 1.5/2.0mm System Instruments Set

98370Monoloc Locking Plates 1.5/2.0mm Instruments Set98371Monoloc Locking Plates 1.5/2.0mm Instruments Set (empty)-PPSU1216100Handle, with Mini Quick Coupling2217490Bending-Cutting Pliers1217500/217510Universal Pliers, pointed for Plate 1.0 to 2.0 (left/right)101123Reduction Forceps, pointed101215Retractor, width 6mm101216Retractor, width 15mm101218Sharp Hook101232Periosteal Elevator, width 5mm101232Periosteal Elevator, width 3mm1217530Screw Holding Forceps for 1.0 to 2.0 Cortex Screws1216180Plate Holding Instruments, 4.01
98371Monoloc Locking Plates 1.5/2.0mm Instruments Set (empty)-PPSU1216100Handle, with Mini Quick Coupling2217490Bending-Cutting Pliers1217500/217510Universal Pliers, pointed for Plate 1.0 to 2.0 (left/right)101123Reduction Forceps, pointed101215Retractor, width 6mm101216Retractor, width 8mm101217Retractor, width 15mm101218Sharp Hook101232Periosteal Elevator, width 3mm101232Periosteal Elevator, width 3mm1217530Screw Holding Forceps for 1.0 to 2.0 Cortex Screws1216180Plate Holding Instruments, 4.01
216100Handle, with Mini Quick Coupling2217490Bending-Cutting Pliers1217500/217510Universal Pliers, pointed for Plate 1.0 to 2.0 (left/right)101123Reduction Forceps, pointed101215Retractor, width 6mm101216Retractor, width 8mm101217Retractor, width 15mm101218Sharp Hook101229Periosteal Elevator, width 3mm101232Periosteal Elevator, width 3mm1217530Screw Holding Forceps for 1.0 to 2.0 Cortex Screws1216180Plate Holding Instruments, 4.01
217490Bending-Cutting Pliers1217500/217510Universal Pliers, pointed for Plate 1.0 to 2.0 (left/right)101123Reduction Forceps, pointed101215Retractor, width 6mm101216Retractor, width 8mm101217Retractor, width 15mm101218Sharp Hook101232Periosteal Elevator, width 3mm101232Periosteal Elevator, width 3mm1217530Screw Holding Forceps for 1.0 to 2.0 Cortex Screws1216180Plate Holding Instruments, 4.01
217500/217510Universal Pliers, pointed for Plate 1.0 to 2.0 (left/right)101123Reduction Forceps, pointed101215Retractor, width 6mm101216Retractor, width 8mm101217Retractor, width 15mm101218Sharp Hook101219Periosteal Elevator, width 5mm101232Periosteal Elevator, width 3mm1217530Screw Holding Forceps for 1.0 to 2.0 Cortex Screws1216180Plate Holding Instruments, 4.01
01123Reduction Forceps, pointed101215Retractor, width 6mm101216Retractor, width 8mm101217Retractor, width 15mm101218Sharp Hook101219Periosteal Elevator, width 5mm101232Periosteal Elevator, width 3mm1217530Screw Holding Forceps for 1.0 to 2.0 Cortex Screws1216180Plate Holding Instruments, 4.01
01215Retractor, width 6mm101216Retractor, width 8mm101217Retractor, width 15mm101218Sharp Hook101219Periosteal Elevator, width 5mm101232Periosteal Elevator, width 3mm1217530Screw Holding Forceps for 1.0 to 2.0 Cortex Screws1216180Plate Holding Instruments, 4.01
01216Retractor, width 8mm101217Retractor, width 15mm101218Sharp Hook101219Periosteal Elevator, width 5mm101232Periosteal Elevator, width 3mm1217530Screw Holding Forceps for 1.0 to 2.0 Cortex Screws1216180Plate Holding Instruments, 4.01
01217Retractor, width 15mm101218Sharp Hook101219Periosteal Elevator, width 5mm101232Periosteal Elevator, width 3mm1217530Screw Holding Forceps for 1.0 to 2.0 Cortex Screws1216180Plate Holding Instruments, 4.01
01218Sharp Hook101219Periosteal Elevator, width 5mm101232Periosteal Elevator, width 3mm1217530Screw Holding Forceps for 1.0 to 2.0 Cortex Screws1216180Plate Holding Instruments, 4.01
01219Periosteal Elevator, width 5mm101232Periosteal Elevator, width 3mm1217530Screw Holding Forceps for 1.0 to 2.0 Cortex Screws1216180Plate Holding Instruments, 4.01
01232Periosteal Elevator, width 3mm1217530Screw Holding Forceps for 1.0 to 2.0 Cortex Screws1216180Plate Holding Instruments, 4.01
217530Screw Holding Forceps for 1.0 to 2.0 Cortex Screws1216180Plate Holding Instruments, 4.01
216180Plate Holding Instruments, 4.01
216190Plate Holding Instruments, 5.01
98567 Holding Forceps for Finger Plates, ratchet lock 1
98086Depth Gauge for Screws Φ 1.5 mm, measuring range up to 24 mm1
98087Depth Gauge for ScrewΦ 2.0 mm, measuring range up to 40 mm1
161210 Cutter 1

1.5mm implants tray

Code	Product Description	Qty
98372	1.5mm implants tray (Empty)	1
98091	Drill Bit Φ 1.1mm, with Quick Coupling	2
98085	Drill Bit Φ 1.5mm, with Quick Coupling	2
98080	Countersink for Screws Φ1.3 and 1.5mm	1
98083	Double Drill Guide 1.5/1.1	1
98094	LCP Drill Sleeve 1.5, for Drill Bits Φ1.1 mm	2
98100	Screwdriver Shaft Stardrive 1.5, T4, self-holding, with Quick Coupling	2
98098	Holding Sleeve 1.5	1
98092	Bending Pin for LCP Plates 1.5, with thread	2

Product No.	Size	Material	Model No
39475520	2x20mm	ТА	SSQ01
39475522	2x22mm	ТА	SSQ01
39475524	2x24mm	ТА	SSQ01
39475526	2x26mm	ТА	SSQ01
39475528	2x28mm	ТА	SSQ01
39475530	2x30mm	ТА	SSQ01
39475532	2x32mm	ТА	SSQ01
39475534	2x34mm	ТА	SSQ01
39475538	2x38mm	TA	SSQ01





Instruments Set

1.5mm implants tray

Code	Product Description	Qty
98372	1.5mm implants tray (Empty)	1
98091	Drill Bit Φ 1.1mm, with Quick Coupling	2
98085	Drill Bit Φ 1.5mm, with Quick Coupling	2
98080	Countersink for Screws Ф1.3 and 1.5mm	1
98083	Double Drill Guide 1.5/1.1	1
98094	LCP Drill Sleeve 1.5, for Drill Bits Φ1.1 mm	2
98100	Screwdriver Shaft Stardrive 1.5, T4, self-holding, with Quick Coupling	2
98098	Holding Sleeve 1.5	1
98092	Bending Pin for LCP Plates 1.5, with thread	2

2.0mm implants tray

Code	Product Description
98373	2.0mm implants tray (Empty)
98085	Drill Bit Φ1.5mm, with Quick Coupling
98082	Drill Bit Φ2.0mm, with Quick Coupling
98093	Countersink 1.5 to 2.4, with Quick Coupling
98095	Double Drill Guide 2.0/1.5 for Drill Bits 2.0 and 1.5mm
98096	LCP Drill Sleeve 2.0, with Scale, for Drill Bits Φ1.5 mm
98099	Screwdriver Shaft Stardrive 2.0, with Quick Coupling
98097	Holding Sleeve 2.0
98081	Bending Pin for LCP Plates 2.0

Please contact your sales representative for final product list.



98370







Instruments





28



98567

98097



