# GRAFTON® and GRAFTON PLUS® M708348B324S Rev. B DeminerIALIZED Bone Matrix (DBM)

GRAFTON AND GRAFTON PLUS® DEMINERALIZED BONE MATRIX (DBM)

GRAFTON® VE GRAFTON® PLUS MINERALLERDEN ARINDIRILMIŞ KEMİK MATRİKSİ (DBM)

**Medtronic Sofamor Danek USA, Inc.**  
1800 Pyramid Place  
Memphis, TN 38132  
Telephone 800 933 2635 (In U.S.A.)  
901 396 3133 (Outside U.S.A.)  
Fax 901 396 0356

**DISTRIBUTED BY:**  
**Spinalgraft Technologies LLC**  
4340 Swinnea Road, Suite 39  
Memphis, TN 38118

2013-04-22

**Medtronic**

ENGLISH

GRAFTON® AND GRAFTON PLUS® DEMINERALIZED BONE MATRIX (DBM)

READ BEFORE USE

THIS GRAFT IS DERIVED FROM HUMAN TISSUE WHICH WAS GENEROUSLY DONATED SO THAT OTHERS MAY BENEFIT.

EACH UNIT IS INTENDED FOR SINGLE PATIENT, SINGLE USE ONLY.

CAUTION: RESTRICTED TO USE BY A PHYSICIAN OR DENTIST.

NO ADDITIONAL STERILIZATION STEPS ARE TO BE PERFORMED BY THE USER.

## DESCRIPTION

GRAFTON® DBM and GRAFTON PLUS® DBM contain demineralized human bone tissue combined with an inert additive to yield a demineralized bone matrix (DBM) allograft product having a particular physical form and/or handling property. GRAFTON® DBM Gel and GRAFTON PLUS® DBM Paste are produced from a powder form of DBM, while GRAFTON® DBM Flex (which includes A-FLEX®), Putty, Matrix, CRUNCH® and Orthoblend are produced from a fiber form of DBM. GRAFTON® DBM CRUNCH also contains demineralized bone chips/cubes, while GRAFTON® DBM Orthoblend contains non-demineralized cancellous bone chips in addition to the DBM. GRAFTON® DBM and GRAFTON PLUS® DBM products are malleable or pliable and can be molded or cut into various sizes and shapes according to the intended implant site.

This DBM allograft product was prepared from human bone tissue recovered from a cadaveric donor using aseptic surgical techniques and microbiologically tested during recovery. The tissue was further processed under aseptic conditions and was treated with antibiotics (gentamicin), cleaned using 70% alcohol, processed with a surfactant, washed with purified water, and sonicated. Subsequent demineralization of the bone tissue (using the D-MIN® proprietary demineralization process) to produce the DBM in this product was performed so that the resulting bone matrix has a calcium phosphate content level that meets current American Association of Tissue Bank (AATB) standards. The demineralized bone matrix (along with the cancellous bone chips for GRAFTON® DBM Orthoblend) was combined with USP anhydrous glycerol (GRAFTON® DBM allografts) or a starch carrier (GRAFTON PLUS® DBM allografts) to form the final allograft product. The final product in packaged form was tested for sterility according to the procedures in the current U.S. Pharmacopoeia and for endotoxins using a qualified test method.

GRAFTON® DBM and GRAFTON PLUS® DBM are demineralized bone allograft products that are osteoconductive as well as osteoinductive in an athymic rat assay.

GRAFTON® DBM and GRAFTON PLUS® DBM are prepared via a proprietary processing method of Medtronic that has been validated to consistently produce DBM that is osteoinductive in an athymic rat assay. Product and process consistency are confirmed via ongoing testing of GRAFTON® DBM and GRAFTON PLUS® DBM finished product for osteoinductivity in this validated athymic rat assay utilizing a five-point linear scale (0,1,2,3,4) to score bone formation at 28 days post implantation\*. This bone forming activity exhibited by GRAFTON® DBM and GRAFTON PLUS® DBM in this athymic rat surrogate assay should not be interpreted as a predictor of clinical performance.

\*Edwards, J.T., PhD, Diegmann, M.H., MS, Scarborough, N.L., PhD.: Osteoinduction of Human Demineralized Bone: Characterization in a Rat Model. *Clinical Orthopaedics*, December, 1998, Vol 357.

GRAFTON® DBM and GRAFTON PLUS® DBM are packaged in ready-to-use form in single patient use containers. The lot number, expiration date, product code, quantity (volume or size), and additional information are listed on the package label.

CAUTION: Federal law (U.S.A.) restricts this product to sale by or on the order of a physician or dentist.

## VIRAL INACTIVATION PROCEDURES

The DBM in GRAFTON® DBM and GRAFTON PLUS® DBM is produced by a proprietary production process that has been validated to inactivate viruses including: HIV-1; hepatitis B virus (duck hepatitis virus as model); hepatitis C virus (bovine diarrhea virus as model), CMV; and Polio virus.† Testing was performed according to the current concepts and study design elements for process validation studies for the removal and/or inactivation of viruses in the production of biopharmaceutical products recommended by the Food and Drug Administration's (FDA) Center for Biologics Evaluation and Research 1, 2, 3 and the European community 4, 5, 6, 7. All studies were performed in certified conformity with Good Laboratory Practice for Nonclinical Laboratory Studies regulations stated in the Code of Federal Regulations (21 CFR § 58).

These viral inactivation procedures were used to further reduce the risk of disease transmission via the use of GRAFTON® DBM and GRAFTON PLUS® DBM allografts beyond the protection provided by donor testing and screening procedures.

The process used to produce the non-demineralized cancellous bone chips in GRAFTON® DBM Orthoblend does not afford the same degree of viral inactivation as the process used to produce the DBM. However, the risk of disease transmission with this tissue component remains low due to multiple safeguards that are rigorously employed, including donor screening, laboratory testing, and material processing.

†Data on file at Medtronic

## INDICATIONS FOR USE

GRAFTON® DBM and GRAFTON PLUS® DBM are intended for use as a bone graft extender, bone graft substitute, and bone void filler in bony voids or gaps of the skeletal system (i.e., spine, pelvis, and extremities) not intrinsic to the stability of the bony structure. The voids or gaps may be surgically created defects or defects created by traumatic injury to the bone.

GRAFTON® DBM (excluding the Orthoblend form) and GRAFTON PLUS® DBM are also intended to be packed into bony voids or gaps to fill and/or augment dental intraseous, oral, and cranio-/maxillofacial defects. These defects may be surgically created osseous defects or osseous defects created from traumatic injury to the bone, including periodontal/infrabony defects; alveolar ridge augmentation (sinusotomy, osteotomy, cystectomy); dental extraction sites (ridge maintenance, implant preparation/ placement); sinus lifts; cystic defects; craniofacial augmentation. GRAFTON® DBM and GRAFTON PLUS® DBM may be used alone in a manner comparable to autogenous bone chips or allograft bone particulate (demineralized freeze dried bone), or they may be mixed with either allograft or autograft bone or bone marrow as a bone graft extender. GRAFTON® DBM and GRAFTON PLUS® DBM are indicated only for bony voids or gaps that are not intrinsic to the stability of the bony structure.

GRAFTON® DBM and GRAFTON PLUS® DBM are absorbed/remodeled and replaced by host bone during the healing process.

Note: The user should consider the fact that GRAFTON® DBM CRUNCH contains demineralized bone chips approximately 3 mm (±1 mm) in determining the appropriateness of this allograft for use in small defects.

## CONTRAINDICATIONS

The following are contraindications for the use of GRAFTON® DBM and GRAFTON PLUS® DBM:

- The presence of infection at the transplantation site.
- Treatment of spinal insufficiency fractures.

## CAUTION

This allograft may contain trace amounts of antibiotics (gentamicin), surfactant, and other processing solutions. Caution should be exercised if the patient is allergic to these antibiotics or chemicals.

GRAFTON PLUS® DBM Paste contains starch. Therefore, caution should be exercised in using GRAFTON PLUS® DBM Paste in a patient with a starch allergy and/or amylase deficiency.

## PRECAUTIONS

Extensive donor blood serum testing, medical and social history screening procedures, and tissue microbiological testing have been used in the qualification of all tissue donors. Despite the viral inactivation and extensive tissue donor selection and qualification processes used in providing this tissue graft, transmission of an infectious disease through the use of this tissue graft is still possible. Bacterial infection at the graft site may also occur. Any adverse outcomes potentially attributable to GRAFTON® DBM or GRAFTON PLUS® DBM must be reported promptly to Medtronic

Adequate fixation should be used to stabilize the implant site during bone formation and healing in bony voids or gaps of the skeletal system (i.e., spine, pelvis, and extremities).

If injecting GRAFTON® DBM or GRAFTON PLUS® DBM into the defect site, precaution should be taken not to:

- over-pressurize the delivery device, as this may lead to extrusion of the device beyond the site of its intended application and damage to the surrounding tissues.
- over-pressurize the defect site, as this may lead to fat embolization or embolization of the device material into the bloodstream.

When used as a bone graft extender in bony voids or gaps of the skeletal system (i.e., spine, pelvis, and extremities), GRAFTON PLUS® DBM Paste is intended for use only with autograft, not other allograft. Recommended ratios of GRAFTON PLUS® DBM Paste to autograft as a bone graft extender are 1:1 or 2:1.

## DONOR SCREENING AND TESTING

Prior to donation, the donor's blood, tissues, and medical/social history were screened for medical conditions or disease processes that would contraindicate the donation of tissues in accordance with current FDA regulations and standards established by the American Association of Tissue Banks. The donor's medical/social history was also screened for HIV, Hepatitis, and CJD/vCJD high risk factors in accordance with current United States Public Health Services Recommendations and FDA Federal Regulations and Guidance Documents.

Testing of donor blood and tissue samples began at the site of recovery and continued into processing. Donor blood samples taken at the time of recovery were tested for communicable disease by a laboratory certified to perform such testing on human specimens under the Clinical Laboratory Improvement Amendments of 1988 (42 U.S.C. 263a) and 42 CFR 493, using FDA licensed tests including:

- HBsAg (Hepatitis B Surface Antigen)
- HBe-IgM/IgG (Hepatitis B Total Core Antibody)
- HCV (Hepatitis C Antibody)
- HIV 1/2-Ab (Antibody to Human Immunodeficiency Virus Types 1 and 2)
- HTLV I/II (Human T-Lymphotropic Virus Types I and II)
- RPR/STS or Equivalent (Syphilis Detection)

In addition to the above tests, donors recovered prior to March 9, 2005, were tested for HIV-1 by p24 antigen testing or DNA Polymerase Chain Reaction (PCR) or Transcription Mediated Amplification (TMA) Nucleic Acid Amplification Testing (NAT). Donors recovered after March 9, 2005, were tested using Human Immunodeficiency Virus Type 1/Nucleic Acid Amplification Testing (HIV1 NAT) and Hepatitis C Virus/Nucleic Acid Amplification Testing (HCV NAT).

The results of all the relevant communicable disease tests referenced above were found to be negative.

The communicable disease test results, together with the informed consent, medical history interview, physical assessment, available medical records (to include previous medical history, laboratory test results, autopsy and coroner reports, if performed), and information obtained from any source or records which may pertain to donor eligibility have been evaluated. Based on this evaluation, the donor was determined to meet donor eligibility criteria current at the time of recovery. The donor eligibility criteria used to screen this donor are in compliance with FDA regulations published in 21 CFR Part 1270 "Human Tissue Intended for Transplantation" and/or Part 1271 "Human Cells, Tissues, and Cellular and Tissue-Based Products", as applicable.

Donor eligibility was determined by one of the following tissue banks:

American Tissue Services Foundation,  
Edmond, OK 73013

Musculoskeletal Transplant Foundation,  
Edison, NJ 08837

Community Tissue Services  
Dayton, OH 45402

LifeNet Health  
Virginia Beach, VA 23453

RTI Biologics  
Alachua, FL 32615

The names and addresses of the testing laboratories, the listing and interpretation of all required communicable disease tests, a listing of the documents reviewed as part of the relevant medical records, and the name of the person or establishment determining the eligibility of this human tissue are on file at Medtronic, Eatontown, NJ and are available upon request.

The final tissue allograft product was released by Medtronic based on the initial donor eligibility determination and on a post-processing review and determination that the product met all processing requirements and specifications.

This tissue allograft product has been released for transplantation.

## STERILITY

GRAFTON® DBM and GRAFTON PLUS® DBM tissue allograft products have been aseptically processed and tested for sterility, as indicated by the package label and as explained below. Do not subject this allograft to additional disinfection or sterilization procedures.

Tissue allografts labeled as "Aseptically Processed, Passes USP Sterility Tests" or "Sterile" were aseptically processed and then tested for sterility according to the procedures in the current U.S. Pharmacopoeia.

Medtronic may use low dose gamma irradiation as an adjunct to aseptic processing to reduce bioburden. GRAFTON® DBM and GRAFTON PLUS® DBM package labels containing "Tissue Gamma Irradiated" indicate that low dose (1.0 – 1.8 megarads) gamma irradiation was used as a means of reducing the bioburden on the donor tissue.

## TISSUE TRACKING

Federal (USA) regulations under 21 CFR 1271 establish requirements for the tracking of human tissue products. In accordance with these regulations, the package label of each GRAFTON® DBM and GRAFTON PLUS® DBM unit distributed by Medtronic bears a lot number that serves as a distinct identification code that is recorded in Medtronic's distribution records for purposes of tracking the tissue to the consignee or user/tissue transplant facility. This lot number should be recorded in the user/tissue transplant facility's records and in the tissue recipient's medical record, along with the following information:

1. Description of Tissue
2. Lot Number (Donor ID)
3. Product Code
4. Expiration Date
5. Quantity Implanted
6. Antibiotics Used
7. Description of Procedure
8. Date and Time of Procedure
9. Surgeon Name
10. Any Other Pertinent Information

If for any reason the tissue allograft is opened and not utilized, it should be disposed of properly or returned to Medtronic. Document the reason for the tissue not being utilized.

For European Economic Area audience only: These tissue tracking records shall be maintained for 30 years after clinical use. In case product traceability to recipient is at risk, the customer shall secure the transfer of the records to another entity (preferably a Tissue Establishment or Organ Bank) in order to secure continued traceability. Medtronic BV shall be informed of such a transfer of records.

## INSTRUCTIONS FOR USE

GRAFTON® DBM and GRAFTON PLUS® DBM have been tested for sterility. Do not subject this product to additional disinfection or sterilization procedures. The contents of an individual GRAFTON® DBM or GRAFTON PLUS® DBM container are intended for single patient use only. Do not use the contents of any container for multiple patients. Empty or partially used containers should be disposed of in accordance with recognized procedures for discarding medical waste materials.

Preparation of the bone graft bed is important for graft incorporation and bone formation, as are other factors such as blood supply, source of marrow elements, loading, stability, and absence of infection at the graft site. The volume of graft material used in each procedure is determined by the judgment of the clinician.

GRAFTON® DBM (Gel, Putty) and GRAFTON PLUS® DBM may be injected into the defect site (see Precautions and Contraindications).

Before Usage: Examine Product Package – Do Not Use This Product If:

- Any of the package materials or contents appear to be missing, tampered with or damaged;
- The package label or identifying bar code is illegible or missing
- The expiration date shown on the package label has passed
- If any of the above conditions exist or are suspected, this allograft should not be used.

NOTE: Once a package seal has been compromised, the tissue should be either transplanted, if appropriate, or otherwise properly discarded.

Opening Instructions

1. Peel open the outer pouch using proper sterile technique.
2. Pass sterile contents into sterile field.
3. For syringe type containers, remove and discard cap prior to use.

### GRAFTON® DBM Gel Preparation for Use

KEEP THE SURGICAL SITE AS DRY AS POSSIBLE. Do not irrigate during or after placement of GRAFTON® DBM Gel. Fluids such as water, saline, or blood may alter the consistency and handling characteristics of GRAFTON® DBM Gel.

If GRAFTON® DBM Gel is used in combination with autologous bone tissues, other forms of allografts or other bone grafting materials, it is recommended to remove from these graft materials any excess fluids (including rehydration solution) before the other grafts come into contact with GRAFTON® DBM Gel.

### GRAFTON® DBM Putty Preparation for Use

GRAFTON® DBM Putty REQUIRES NO REHYDRATION prior to use. If desired, a small amount of fluid such as blood, sterile water, or sterile saline may be added to GRAFTON® DBM Putty in order to adjust its consistency or handling characteristics. Also, simply kneading the GRAFTON® DBM Putty may enhance the pliability and cohesiveness of the product.

### GRAFTON® DBM Flex (including A-FLEX) Preparation for Use

GRAFTON® DBM Flex may BE REHYDRATED prior to use to attain porosity and full flexibility. Place in sterile saline or blood until desired consistency is achieved.

### GRAFTON® DBM CRUNCH and Orthoblend Preparation for Use

GRAFTON® DBM CRUNCH and Orthoblend REQUIRE NO REHYDRATION prior to use. If desired, a small amount of fluid such as bone marrow aspirate, blood, sterile water, or sterile saline may be added to GRAFTON® DBM CRUNCH or Orthoblend in order to adjust the consistency or handling characteristics. Also, simply kneading the GRAFTON® DBM CRUNCH or Orthoblend may enhance the pliability and cohesiveness of the allograft.

### GRAFTON® DBM Matrix Preparation for Use

GRAFTON® DBM Matrix MAY BE REHYDRATED prior to use. If desired, a small amount of fluid such as blood, sterile water, or sterile saline may be added to GRAFTON® DBM Matrix in order to adjust its consistency and handling characteristics.

### GRAFTON PLUS® DBM Paste Preparation for Use

GRAFTON PLUS® DBM Paste REQUIRES NO REHYDRATION prior to use. If desired, a small amount of fluid such as blood, sterile water, or sterile saline may be added to GRAFTON PLUS® DBM Paste in order to adjust its consistency or handling characteristics. Also, simply kneading the GRAFTON PLUS® DBM Paste may enhance the pliability and cohesiveness of the product.

Note: GRAFTON PLUS® DBM Paste will turn purple or black in color upon contact with iodine due to the presence of starch in the product. This color change merely indicates that the product has absorbed some iodine and does not present any safety issues beyond those associated with the iodine itself. If contact of the product with iodine should occur, the surgeon should use discretion regarding use of the product just as with any other grafting material that would come into contact with iodine.

## STORAGE

Refer to product package label for specific storage conditions. It is the responsibility of the transplant facility or clinician to maintain the tissue intended for transplantation in the appropriate recommended storage conditions prior to transplant.

## RETURNS

For any product returns to Medtronic, a Return Authorization Number is required prior to returning this product. Refer to Medtronic's Return Policy.

For European Economic Area audience only: The customer shall notify serious adverse events (meaning any untoward occurrence associated with the procurement, testing, processing, storage, and distribution of tissues and cells that might lead to the transmission of a communicable disease, to death or life-threatening, disabling or incapacitating conditions for patients or which might result in, or prolong, hospitalization or morbidity) or serious adverse reactions (meaning an unintended response, including a communicable disease, in the donor or in the recipient associated with the procurement or human application of tissues and cells that is fatal, life-threatening, disabling, incapacitating or which results in, or prolongs, hospitalization or morbidity) within 1 day to their local Medtronic representative.

## MRI INFORMATION

GRAFTON® DBM and GRAFTON PLUS® DBM are MR Safe.

GRAFTON® DBM and GRAFTON PLUS® DBM are nonconducting or nonmagnetic items which pose no known hazards in all MR environments for magnetically induced displacement force and magnetically induced torque. In addition, GRAFTON® DBM and GRAFTON PLUS® DBM are not susceptible to heating due to RF (radio frequency) fields. As such, GRAFTON® DBM and GRAFTON PLUS® DBM can justifiably be labeled as MR-Safe per ASTM F2503.

## REFERENCES

Standards for Tissue Banking (current version), American Association of Tissue Banks, Arlington, VA.

Current Policies and Procedures of Medtronic, Eatontown, N.J.

21CFR1270 titled "Human Tissue Intended for transplantation"

21CFR1271 titled "Human Cells, Tissues, and Cellular and Tissue-Based Products"

PHS Guidelines for Preventing Transmission of HIV through Transplantation of Human Tissue and Organs, MMWR 1994:43, 1-17.

PHS Guideline for Screening Donors of Blood, Plasma, Organs, Tissue and Semen for Evidence of Hepatitis B and Hepatitis C, MMWR 1991:40, 1-17.

FDA Recommendations to Blood Establishments for "Deferral of Current and Recent Inmates of Correctional Institutions as Donors of Whole Blood, Blood Components, Source Leukocytes, and Source Plasma," 6/8/95.

FDA Recommendations to Blood Establishments for "Precautionary Measures to Further Reduce the Possible Risk of Transmission of Creutzfeldt-Jakob Disease By Blood and Blood products," 8/8/95.

Notes

1. Center for Biologics Evaluation and Research, "Points to Consider in the Characterization of Cell Lines Used to Produce Biologicals" (Food and Drug Administration) 1993.
2. Center for Biologics Evaluation and Research, "Points to Consider in the Manufacture and Testing of Monoclonal Antibody Products for Human Use" (Food and Drug Administration) 1994.
3. Kozak RW. Viral Removal and Inactivation Issues for Biological Products. Proceedings of the 1991 Technical Program Pharmaceutical and Cosmetic Industries Exposition and Conference. 1991: 253-260.
4. Committee for Proprietary Medicinal Products: Ad Hoc Working Party on Biotechnology/Pharmacy and Working Party on Safety Medicines, Note for Guidance: "Validation of Virus Removal and Inactivation Procedures", Biologicals 1991; 19-247-251.
5. Committee for Proprietary Medicinal Products: EEC Council Directive 89/381: "Medicinal Products Derived from Human Plasma" (Revised Draft 1995).
6. Committee for Proprietary Medicinal Products: 1995 Revised CPMP Guidelines. Virus Validation Studies: The design, contribution, and interpretation of studies validating the inactivation and removal of viruses (revised).
7. ICH Viral Safety Document Draft April 1995: Viral Safety Evaluation of Biotechnology Products Derived from Cell Lines of Human or Animal Origin.

GRAFTON®, GRAFTON PLUS®, D-MIN®, CRUNCH®, A-FLEX®, ADVANCING OSTEOBIOLOGIC SCIENCES® and MEDTRONIC® are registered trademarks of Medtronic, Inc.

Covered by one or more of U.S. Patent Nos. 5,314,476; 5,484,601; 5,507,813; 5,607,269; 5,676,146 C2; and 7,163,691, and foreign patents. Other U.S. and foreign patents pending.

©2013 Medtronic Sofamor Danek USA, Inc. All rights reserved.

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