


VARIOLOC

Locking Compression Plate System • Catalog 2016



THERAPY OPTIMIZED
EXPECT EXTRAORDINARY

TABLE OF CONTENT


Introduction	VARILOC Locking Compression Plate System	2
	FEATURES & BENEFITS	2
Know Your Options by Anatomy	Shoulder	3
	Hand & Wrist	7
	Hip/Pelvis	14
	Knee	15
	Foot & Ankle	17
	Bone Shaft	18
	Screws	21
Instrument	Instruments	32

Warning
Please contact your sales representative for detailed product list.

INTRODUCTION

VARILOC Locking Compression Plate system inherited the advantages of locking plate which provides angle stability to achieve better fixation in metaphyseal fractures and osteoporosis. Thanks to its locking mechanism, it offers better preservation of blood supply comparing to conventional plate and screw system. The combination hole design in the plate shaft combines function of dynamic compression hole and threaded locking hole which provides the flexibility of cortex screw or locking screw fixation. As innovation in this system, Variable Angle Locking Screw and Cannulated Locking Screws offer surgeons with more accurate screw insertion angle and optimal intercondylar compression.

FEATURES & BENEFITS



KNOW YOUR OPTIONS BY ANATOMY

Shoulder

Clavicular S Reconstruction Locking Plate 3.5



Specification	
Thickness (mm)	2.8
Width (mm)	10
Hole Space (mm)	12
Screws	3.5 cortex screw
	3.5 locking screw
	3.5 VA screw

Code	Direction	Size	Length (mm)	Model No.	Material
24256106	R	6H	61.8	SDXZ27	PT
24256107	R	7H	77.1	SDXZ27	PT
24256108	R	8H	91.3	SDXZ27	PT
24256109	R	9H	92	SDXZ27	PT
24256110	R	10H	103.9	SDXZ27	PT
24256112	R	12H	127.7	SDXZ27	PT
24256206	L	6H	61.8	SDXZ28	PT
24256207	L	7H	77.1	SDXZ28	PT
24256208	L	8H	91.3	SDXZ28	PT
24256209	L	9H	92	SDXZ28	PT
24256210	L	10H	103.9	SDXZ28	PT
24256212	L	12H	127.7	SDXZ28	PT



KNOW YOUR OPTIONS BY ANATOMY

Shoulder

Reconstruction Clavicular Locking Compression Plate 3.5



Specification	
Thickness (mm)	3.2
Width (mm)	11
Screws	3.5 cortex screw
	3.5 locking screw

Code	Direction	Size	Length (mm)	Model No.	Material
24257106	R	6H	103.2	SDXG22	PT
24257107	R	7H	113	SDXG22	PT
24257108	R	8H	122.5	SDXG22	PT
24257206	L	6H	103.2	SDXG23	PT
24257207	L	7H	113	SDXG23	PT
24257208	L	8H	122.5	SDXG23	PT

Anterolateral Clavicular Locking Compression Plate Plate, 2.7/3.5



Specification	
Thickness (mm)	3.3
Width (mm)	10.5
Head Screws	2.7 locking screw
	2.7 VA locking screw
Shaft Screws	3.5 cortex screw
	3.5 locking screw

Code	Direction	Size	Length (mm)	Model No.	Material
24258110	R	7+3H	68	SDXG07	PT
24258111	R	7+4H	80	SDXG07	PT
24258112	R	7+5H	92	SDXG07	PT
24258113	R	7+6H	104	SDXG07	PT
24258114	R	7+7H	116	SDXG07	PT
24258115	R	7+8H	128	SDXG07	PT
24258210	L	7+3H	68	SDXG08	PT
24258211	L	7+4H	80	SDXG08	PT
24258212	L	7+5H	92	SDXG08	PT
24258213	L	7+6H	104	SDXG08	PT
24258214	L	7+7H	116	SDXG08	PT
24258215	L	7+8H	128	SDXG08	PT

KNOW YOUR OPTIONS BY ANATOMY

Shoulder

Clavicle Hook Locking Compression Plate II, 12mm

					
Specification					
Thickness (mm)		3			
Width (mm)		11			
Head Screws		2.7 VA locking screw 2.7 locking screw			
Shaft Screws		3.5 locking screw 3.5 cortex screw			
Code	Direction	Size	Length (mm)	Model No.	Material
24282102	R	3+2H	44	SDXG26	PT
24282103	R	3+3H	54	SDXG26	PT
24282104	R	3+4H	64	SDXG26	PT
24282105	R	3+5H	74	SDXG26	PT
24282106	R	3+6H	84	SDXG26	PT
24282107	R	3+7H	94	SDXG26	PT
24282202	L	3+2H	44	SDXG25	PT
24282203	L	3+3H	54	SDXG25	PT
24282204	L	3+4H	64	SDXG25	PT
24282205	L	3+5H	74	SDXG25	PT
24282206	L	3+6H	84	SDXG25	PT
24282207	L	3+7H	94	SDXG25	PT

Clavicle Hook Locking Compression Plate II, 15mm

					
Specification					
Thickness (mm)		3			
Width (mm)		11			
Head Screws		2.7 VA locking screw 2.7 locking screw			
Shaft Screws		3.5 locking screw 3.5 cortex screw			
Code	Direction	Size	Length (mm)	Model No.	Material
24281102	R	3+2H	45	SDXG26	PT
24281103	R	3+3H	55	SDXG26	PT
24281104	R	3+4H	65	SDXG26	PT
24281105	R	3+5H	75	SDXG26	PT
24281106	R	3+6H	85	SDXG26	PT
24281107	R	3+7H	95	SDXG26	PT
24281202	L	3+2H	45	SDXG25	PT
24281203	L	3+3H	55	SDXG25	PT
24281204	L	3+4H	65	SDXG25	PT
24281205	L	3+5H	75	SDXG25	PT
24281206	L	3+6H	85	SDXG25	PT
24281207	L	3+7H	95	SDXG25	PT

KNOW YOUR OPTIONS BY ANATOMY

Shoulder

Clavicle Hook Locking Compression Plate II, 18mm

					
Specification					
Thickness (mm)		3			
Width (mm)		11			
Head Screws		2.7 VA locking screw 2.7 locking screw			
Shaft Screws		3.5 locking screw 3.5 cortex screw			
Code	Direction	Size	Length (mm)	Model No.	Material
24283102	R	3+2H	46	SDXG26	PT
24283103	R	3+3H	56	SDXG26	PT
24283104	R	3+4H	66	SDXG26	PT
24283105	R	3+5H	76	SDXG26	PT
24283106	R	3+6H	86	SDXG26	PT
24283107	R	3+7H	96	SDXG26	PT
24283202	L	3+2H	46	SDXG25	PT
24283203	L	3+3H	56	SDXG25	PT
24283204	L	3+4H	66	SDXG25	PT
24283205	L	3+5H	76	SDXG25	PT
24283206	L	3+6H	86	SDXG25	PT
24283207	L	3+7H	96	SDXG25	PT

KNOW YOUR OPTIONS BY ANATOMY


Hand & Wrist

Distal Radius Two-Column Volar Plate I, 2.4/2.7



Specification	
Thickness (mm)	1.7
Head Width (mm)	19.5
Head Screws	2.4 locking screw 2.7 locking screw 2.4 VA locking screw 2.7 VA locking screw
Shaft Screws	2.4 cortex screw 2.7 cortex screw 2.7 locking screw

Code	Direction	Size	Length (mm)	Model No.	Material
24278102	R	7+2H	54	SDDG16	PT
24278103	R	7+3H	63	SDDG16	PT
24278104	R	7+4H	72	SDDG16	PT
24278105	R	7+5H	81	SDDG16	PT
24278202	L	7+2H	54	SDDG15	PT
24278203	L	7+3H	63	SDDG15	PT
24278204	L	7+4H	72	SDDG15	PT
24278205	L	7+5H	81	SDDG15	PT



KNOW YOUR OPTIONS BY ANATOMY

Hand & Wrist

Distal Radius Two-Column Volar Plate II, 2.4/2.7



Specification	
Thickness (mm)	2
Head Width (mm)	22
Head Screws	2.4 locking screw 2.7 locking screw 2.4 VA locking screw 2.7 VA locking screw
Shaft Screws	2.4 cortex screw 2.7 cortex screw 2.7 locking screw

Code	Direction	Size	Length (mm)	Model No.	Material
24279102	R	7+2H	54	SDDG16	PT
24279103	R	7+3H	63	SDDG16	PT
24279104	R	7+4H	72	SDDG16	PT
24279105	R	7+5H	81	SDDG16	PT
24279202	L	7+2H	54	SDDG15	PT
24279203	L	7+3H	63	SDDG15	PT
24279204	L	7+4H	72	SDDG15	PT
24279205	L	7+5H	81	SDDG15	PT

Distal Radius Two-Column Volar Plate III, 2.4/2.7



Specification	
Thickness (mm)	2.3
Head Width (mm)	25.5
Head Screws	2.4 locking screw 2.7 locking screw 2.4 VA locking screw 2.7 VA locking screw
Shaft Screws	2.4 cortex screw 2.7 cortex screw 2.7 locking screw

Code	Direction	Size	Length (mm)	Model No.	Material
24280102	R	7+2H	54	SDDG16	PT
24280103	R	7+3H	63	SDDG16	PT
24280104	R	7+4H	72	SDDG16	PT
24280105	R	7+5H	81	SDDG16	PT
24280202	L	7+2H	54	SDDG15	PT
24280203	L	7+3H	63	SDDG15	PT
24280204	L	7+4H	72	SDDG15	PT
24280205	L	7+5H	81	SDDG15	PT

KNOW YOUR OPTIONS BY ANATOMY

Hand & Wrist

Distal Radius Volar Column Locking Compression Plate II, 2.4/2.7



Specification	
Thickness (mm)	1.8
Width (mm)	8.5
Hole Space (mm)	9
Head Screws	2.4 locking screw 2.7 locking screw 2.4 VA locking screw 2.7 VA locking screw
Shaft Screws	2.4 cortex screw 2.7 cortex screw 2.7 locking screw

Code	Direction	Size	Length (mm)	Model No.	Material
24259103	R	8+3H	46.5	SDXG59	PT
24259104	R	8+4H	55.5	SDXG59	PT
24259105	R	8+5H	64.5	SDXG59	PT
24259106	R	8+6H	73.5	SDXG59	PT
24259107	R	8+7H	82.5	SDXG59	PT
24259108	R	8+8H	91.5	SDXG59	PT
24259109	R	8+9H	100.5	SDXG59	PT
24259110	R	8+10H	109.5	SDXG59	PT
24259111	R	8+11H	118.5	SDXG59	PT
24259112	R	8+12H	127.5	SDXG59	PT
24259203	L	8+3H	46.5	SDXG59	PT
24259204	L	8+4H	55.5	SDXG59	PT
24259205	L	8+5H	64.5	SDXG59	PT
24259206	L	8+6H	73.5	SDXG59	PT
24259207	L	8+7H	82.5	SDXG59	PT
24259208	L	8+8H	91.5	SDXG59	PT
24259209	L	8+9H	100.5	SDXG59	PT
24259210	L	8+10H	109.5	SDXG59	PT
24259211	L	8+11H	118.5	SDXG59	PT
24259212	L	8+12H	127.5	SDXG59	PT

KNOW YOUR OPTIONS BY ANATOMY

Hand & Wrist

Distal Radius Volar Locking Compression Plate, oblique T , 2.4/2.7



Specification	
Thickness (mm)	2
Width (mm)	7.5
Hole Space (mm)	8
Head Screws	2.4 locking screw 2.7 locking screw 2.4 VA locking screw 2.7 VA locking screw
Shaft Screws	2.4 cortex screw 2.7 cortex screw 2.7 locking screw

Code	Direction	Size	Length (mm)	Model No.	Material
24284103	R	5+3H	45	SDXG54	PT
24284104	R	5+4H	53	SDXG54	PT
24284105	R	5+5H	61	SDXG54	PT
24284106	R	5+6H	69	SDXG54	PT
24284107	R	5+7H	77	SDXG54	PT
24284108	R	5+8H	85	SDXG54	PT
24284109	R	5+9H	93	SDXG54	PT
24284110	R	5+10H	101	SDXG54	PT
24284111	R	5+11H	109	SDXG54	PT
24284112	R	5+12H	117	SDXG54	PT
24284203	L	5+3H	45	SDXG53	PT
24284204	L	5+4H	53	SDXG53	PT
24284205	L	5+5H	61	SDXG53	PT
24284206	L	5+6H	69	SDXG53	PT
24284207	L	5+7H	77	SDXG53	PT
24284208	L	5+8H	85	SDXG53	PT
24284209	L	5+9H	93	SDXG53	PT
24284210	L	5+10H	101	SDXG53	PT
24284211	L	5+11H	109	SDXG53	PT
24284212	L	5+12H	117	SDXG53	PT

KNOW YOUR OPTIONS BY ANATOMY

Hand & Wrist

Distal Dorsal Radius Locking Compression Plate, oblique L , 2.4/2.7

					
Specification					
Thickness (mm)		1.8			
Width (mm)		6.5			
Hole Space (mm)		10			
Head Screws		2.4 locking screw 2.7 locking screw 2.4 VA locking screw 2.7 VA locking screw			
Shaft Screws		2.4 cortex screw 2.7 cortex screw 2.7 locking screw			
Code	Direction	Size	Length (mm)	Model No.	Material
24275103	R	3+3H	43	SDXG52	PT
24275104	R	3+4H	52	SDXG52	PT
24275203	L	3+3H	43	SDXG51	PT
24275204	L	3+4H	52	SDXG51	PT

Distal Dorsal Radius T-Locking Compression Plate, 2.4/2.7

				
Specification				
Thickness (mm)		1.8		
Width (mm)		6.5		
Hole Space (mm)		12		
Head Screws		2.4 locking screw 2.7 locking screw 2.4 VA locking screw 2.7 VA locking screw		
Shaft Screws		2.4 cortex screw 2.7 cortex screw 2.7 locking screw		
Code	Size	Length (mm)	Model No.	Material
24276003	3+3H	42	SDXG43	PT
24276004	3+4H	54	SDXG43	PT
24276005	3+5H	66	SDXG43	PT
24276006	3+6H	78	SDXG43	PT
24276007	3+7H	90	SDXG43	PT
24276008	3+8H	102	SDXG43	PT

KNOW YOUR OPTIONS BY ANATOMY

Hand & Wrist

Distal Dorsal Radius Locking Compression Plate II, 2.4/2.7

					
Specification					
Thickness (mm)		1.5			
Width (mm)		6.5			
Hole Space (mm)		9			
Head Screws		2.4 locking screw 2.7 locking screw 2.4 VA locking screw 2.7 VA locking screw			
Shaft Screws		2.4 cortex screw 2.7 cortex screw 2.7 locking screw			
Code	Direction	Size	Length (mm)	Model No.	Material
24277105	R	2+3H	40	SDXG34	PT
24277106	R	2+4H	49	SDXG34	PT
24277205	L	2+3H	40	SDXG35	PT
24277206	L	2+4H	49	SDXG35	PT

Distal Lateral Radius Locking Compression Plate, 2.4/2.7

				
Specification				
Thickness (mm)		2		
Width (mm)		7		
Hole Space (mm)		9		
Head Screws		2.4 locking screw 2.7 locking screw 2.4 VA locking screw 2.7 VA locking screw		
Shaft Screws		2.4 cortex screw 2.7 cortex screw 2.7 locking screw		
Code	Size	Length (mm)	Model No.	Material
24286005	5H	48	SDXG33	PT
24286006	6H	57	SDXG33	PT

KNOW YOUR OPTIONS BY ANATOMY

Hand & Wrist

Distal Lateral Ulna Locking Compression Plate, 2.4/2.7



Specification	
Thickness (mm)	2.3
Width (mm)	6.5
Hole Space (mm)	9
Head Screws	2.4 locking screw 2.7 locking screw 2.4 VA locking screw 2.7 VA locking screw
Shaft Screws	2.4 cortex screw 2.7 cortex screw 2.7 locking screw

Code	Size	Length (mm)	Model No.	Material
24289006	6H	57	SDXG33	PT
24289007	7H	66	SDXG33	PT

KNOW YOUR OPTIONS BY ANATOMY

Hip/Pelvis

Proximal Lateral Femoral Plate 7.0



Specification	
Thickness (mm)	6
Width (mm)	12
Hole Space (mm)	16
Head Screws	7.0 cannulated locking screw
Shaft Screws	5.0 locking screw

Code	Direction	Size	Length (mm)	Model No.	Material
24252101	R	3+1H	59	SDDZ11	PT
24252102	R	3+2H	75	SDDZ11	PT
24252103	R	3+3H	91	SDDZ11	PT
24252201	L	3+1H	59	SDDZ10	PT
24252202	L	3+2H	75	SDDZ10	PT
24252203	L	3+3H	91	SDDZ10	PT



KNOW YOUR OPTIONS BY ANATOMY

Knee

Distal Lateral Femoral Locking Compression Plate 5.0



Specification	
Thickness (mm)	5.8
Width (mm)	16
Hole Space (mm)	20
Head Screws	5.0 locking screw
	5.0 VA locking screw
	5.5 cannulated locking screw
Shaft Screws	4.5 cortex screw
	5.0 locking screw

Code	Direction	Size	Length (mm)	Model No.	Material
24260105	R	7+5H	156	SDDG53	PT
24260107	R	7+7H	196	SDDG53	PT
24260109	R	7+9H	236	SDDG53	PT
24260111	R	7+11H	276	SDDG53	PT
24260113	R	7+13H	316	SDDG53	PT
24260205	L	7+5H	156	SDDG52	PT
24260207	L	7+7H	196	SDDG52	PT
24260209	L	7+9H	236	SDDG52	PT
24260211	L	7+11H	276	SDDG52	PT
24260213	L	7+13H	316	SDDG52	PT



KNOW YOUR OPTIONS BY ANATOMY

Knee

Proximal Lateral Tibia L Locking Compression Plate 3.5



Specification	
Thickness (mm)	4
Width (mm)	14
Hole Space (mm)	20
Head Screws	3.5 locking screw
	3.5 VA locking screw
	4.0 cannulated locking screw
Shaft Screws	3.5 cortex screw
	3.5 locking screw

Code	Direction	Size	Length (mm)	Model No.	Material
24269109	R	4+9H	120	SDDG48	PT
24269110	R	4+10H	140	SDDG48	PT
24269111	R	4+11H	160	SDDG48	PT
24269112	R	4+12H	180	SDDG48	PT
24269113	R	4+13H	200	SDDG48	PT
24269114	R	4+14H	220	SDDG48	PT
24269115	R	4+15H	240	SDDG48	PT
24269209	L	4+9H	120	SDDG49	PT
24269210	L	4+10H	140	SDDG49	PT
24269211	L	4+11H	160	SDDG49	PT
24269212	L	4+12H	180	SDDG49	PT
24269213	L	4+13H	200	SDDG49	PT
24269214	L	4+14H	220	SDDG49	PT
24269215	L	4+15H	240	SDDG49	PT

KNOW YOUR OPTIONS BY ANATOMY

Foot & Ankle

Calcaneus Locking Plate II 3.5



Specification	
Thickness (mm)	1.5
Width (mm)	36
Screws	3.5 locking screw

Code	Size	Length (mm)	Model No.	Material
24263013	13H	60	SDXZ66	PT
24263014	14H	64	SDXZ66	PT
24263015	15H	70	SDXZ66	PT
24263016	16H	78	SDXZ66	PT

Calcaneus Locking Plate III 3.5



Specification	
Thickness (mm)	1.5
Width (mm)	36
Screws	3.5 locking screw

Code	Size	Length (mm)	Model No.	Material
24264060	12H	60	SDXZ67	PT
24264068	12H	68	SDXZ67	PT



KNOW YOUR OPTIONS BY ANATOMY

Bone Shaft

Straight Locking Compression Plate 3.5, thickness 2.0



Specification	
Thickness (mm)	2
Width (mm)	9.5
Hole Space (mm)	13
Head Screws	3.5 cortex screw 3.5 locking screw

Code	Size	Length (mm)	Model No.	Material
24267004	4H	60	SDXG11	PT
24267005	5H	73	SDXG11	PT
24267006	6H	86	SDXG11	PT
24267007	7H	99	SDXG11	PT
24267008	8H	112	SDXG11	PT
24267010	10H	138	SDXG11	PT
24267012	12H	164	SDXG11	PT

Straight Locking Compression Plate 3.5, thickness 3.2



Specification	
Thickness (mm)	3.2
Width (mm)	10.5
Hole Space (mm)	13
Head Screws	3.5 cortex screw 3.5 locking screw

Code	Size	Length (mm)	Model No.	Material
24266004	4H	60	SDXG11	PT
24266005	5H	73	SDXG11	PT
24266006	6H	86	SDXG11	PT
24266007	7H	99	SDXG11	PT
24266008	8H	112	SDXG11	PT
24266009	9H	125	SDXG11	PT
24266010	10H	138	SDXG11	PT
24266011	11H	151	SDXG11	PT
24266012	12H	164	SDXG11	PT
24266014	14H	190	SDXG11	PT
24266016	16H	216	SDXG11	PT

KNOW YOUR OPTIONS BY ANATOMY

Bone Shaft

Straight Locking Compression Plate 3.5, thickness 3.8

		Specification		
		Thickness (mm)	3.8	
		Width (mm)	12	
		Hole Space (mm)	13	
		Head Screws	3.5 cortex screw 3.5 locking screw	
Code	Size	Length (mm)	Model No.	Material
24268004	4H	51	SDXG37	PT
24268005	5H	64	SDXG37	PT
24268006	6H	77	SDXG37	PT
24268007	7H	90	SDXG37	PT
24268008	8H	103	SDXG37	PT
24268009	9H	116	SDXG37	PT
24268010	10H	129	SDXG37	PT
24268011	11H	142	SDXG37	PT
24268012	12H	155	SDXG37	PT
24268014	14H	181	SDXG37	PT
24268016	16H	207	SDXG37	PT

Straight Locking Compression Plate 3.5, thickness 4.0

		Specification		
		Thickness (mm)	4	
		Width (mm)	14	
		Hole Space (mm)	13	
		Head Screws	3.5 cortex screw 3.5 locking screw	
Code	Size	Length (mm)	Model No.	Material
24274004	4H	51	SDXG37	PT
24274005	5H	64	SDXG37	PT
24274006	6H	77	SDXG37	PT
24274007	7H	90	SDXG37	PT
24274008	8H	103	SDXG37	PT
24274009	9H	116	SDXG37	PT
24274010	10H	129	SDXG37	PT
24274011	11H	142	SDXG37	PT
24274012	12H	155	SDXG37	PT
24274014	14H	181	SDXG37	PT
24274016	16H	207	SDXG37	PT

KNOW YOUR OPTIONS BY ANATOMY

Bone Shaft

Straight Locking Compression Plate 4.5/5.0

		Specification		
		Thickness (mm)	5.8	
		Width (mm)	15.5	
		Hole Space (mm)	18	
		Head Screws	4.5 cortex screw 5.0 locking screw	
Code	Size	Length (mm)	Model No.	Material
24265004	4H	71	SDDG39	PT
24265005	5H	89	SDDG39	PT
24265006	6H	107	SDDG39	PT
24265007	7H	125	SDDG39	PT
24265008	8H	143	SDDG39	PT
24265009	9H	161	SDDG39	PT
24265010	10H	179	SDDG39	PT
24265011	11H	197	SDDG39	PT
24265012	12H	215	SDDG39	PT
24265013	13H	233	SDDG39	PT
24265014	14H	251	SDDG39	PT
24265015	15H	269	SDDG39	PT
24265016	16H	287	SDDG39	PT
24265017	17H	305	SDDG39	PT
24265018	18H	323	SDDG39	PT

SCREWS

Cortex Screw, Fully Threaded,Self-tapping, 2.4



Specification	
Thread dia.	Φ2.4 mm
Shaft dia.	Φ1.8 mm
Drill bit	Φ1.8 mm

Code	Size	Model No.	Material
30374006	Φ2.4×6	HAQ06	TA
30374007	Φ2.4×7	HAQ06	TA
30374008	Φ2.4×8	HAQ06	TA
30374009	Φ2.4×9	HAQ06	TA
30374010	Φ2.4×10	HAQ06	TA
30374011	Φ2.4×11	HAQ06	TA
30374012	Φ2.4×12	HAQ06	TA
30374013	Φ2.4×13	HAQ06	TA
30374014	Φ2.4×14	HAQ06	TA
30374016	Φ2.4×16	HAQ06	TA
30374018	Φ2.4×18	HAQ06	TA

Code	Size	Model No.	Material
30374020	Φ2.4×20	HAQ06	TA
30374022	Φ2.4×22	HAQ06	TA
30374024	Φ2.4×24	HAQ06	TA
30374026	Φ2.4×26	HAQ06	TA
30374028	Φ2.4×28	HAQ06	TA
30374030	Φ2.4×30	HAQ06	TA
30374032	Φ2.4×32	HAQ06	TA
30374034	Φ2.4×34	HAQ06	TA
30374036	Φ2.4×36	HAQ06	TA
30374038	Φ2.4×38	HAQ06	TA
30374040	Φ2.4×40	HAQ06	TA

Cortex Screw, Fully Threaded,Self-tapping, 2.7



Specification	
Thread dia.	Φ2.7 mm
Shaft dia.	Φ1.9 mm
Drill bit	Φ2.0 mm

Code	Size	Model No.	Material
30375006	Φ2.7×6	HAQ06	TA
30375007	Φ2.7×7	HAQ06	TA
30375008	Φ2.7×8	HAQ06	TA
30375009	Φ2.7×9	HAQ06	TA
30375010	Φ2.7×10	HAQ06	TA
30375011	Φ2.7×11	HAQ06	TA
30375012	Φ2.7×12	HAQ06	TA
30375013	Φ2.7×13	HAQ06	TA
30375014	Φ2.7×14	HAQ06	TA
30375016	Φ2.7×16	HAQ06	TA
30375018	Φ2.7×18	HAQ06	TA

Code	Size	Model No.	Material
30375020	Φ2.7×20	HAQ06	TA
30375022	Φ2.7×22	HAQ06	TA
30375024	Φ2.7×24	HAQ06	TA
30375026	Φ2.7×26	HAQ06	TA
30375028	Φ2.7×28	HAQ06	TA
30375030	Φ2.7×30	HAQ06	TA
30375032	Φ2.7×32	HAQ06	TA
30375034	Φ2.7×34	HAQ06	TA
30375036	Φ2.7×36	HAQ06	TA
30375038	Φ2.7×38	HAQ06	TA
30375040	Φ2.7×40	HAQ06	TA

SCREWS

Cortex Screw, Fully Threaded,Self-tapping, 3.5



Specification	
Thread dia.	Φ3.5 mm
Shaft dia.	Φ2.7 mm
Drill bit	Φ2.5 mm

Code	Size	Model No.	Material
30387010	Φ3.5×10	HAQ06	TA
30387012	Φ3.5×12	HAQ06	TA
30387014	Φ3.5×14	HAQ06	TA
30387016	Φ3.5×16	HAQ06	TA
30387018	Φ3.5×18	HAQ06	TA
30387020	Φ3.5×20	HAQ06	TA
30387022	Φ3.5×22	HAQ06	TA
30387024	Φ3.5×24	HAQ06	TA
30387026	Φ3.5×26	HAQ06	TA
30387028	Φ3.5×28	HAQ06	TA

Code	Size	Model No.	Material
30387030	Φ3.5×30	HAQ06	TA
30387032	Φ3.5×32	HAQ06	TA
30387034	Φ3.5×34	HAQ06	TA
30387036	Φ3.5×36	HAQ06	TA
30387038	Φ3.5×38	HAQ06	TA
30387040	Φ3.5×40	HAQ06	TA
30387045	Φ3.5×45	HAQ06	TA
30387050	Φ3.5×50	HAQ06	TA
30387055	Φ3.5×55	HAQ06	TA
30387060	Φ3.5×60	HAQ06	TA

Cortex Screw, Fully Threaded,Self-tapping, 4.5



Specification	
Thread dia.	Φ4.5 mm
Shaft dia.	Φ3.0 mm
Drill bit	Φ3.2 mm

Code	Size	Model No.	Material
30386012	Φ4.5×12	HAQ06	TA
30386014	Φ4.5×14	HAQ06	TA
30386016	Φ4.5×16	HAQ06	TA
30386018	Φ4.5×18	HAQ06	TA
30386020	Φ4.5×20	HAQ06	TA
30386022	Φ4.5×22	HAQ06	TA
30386024	Φ4.5×24	HAQ06	TA
30386026	Φ4.5×26	HAQ06	TA
30386028	Φ4.5×28	HAQ06	TA
30386030	Φ4.5×30	HAQ06	TA
30386032	Φ4.5×32	HAQ06	TA
30386034	Φ4.5×34	HAQ06	TA
30386036	Φ4.5×36	HAQ06	TA
30386038	Φ4.5×38	HAQ06	TA
30386040	Φ4.5×40	HAQ06	TA
30386042	Φ4.5×42	HAQ06	TA

Code	Size	Model No.	Material
30386044	Φ4.5×44	HAQ06	TA
30386046	Φ4.5×46	HAQ06	TA
30386048	Φ4.5×48	HAQ06	TA
30386050	Φ4.5×50	HAQ06	TA
30386052	Φ4.5×52	HAQ06	TA
30386054	Φ4.5×54	HAQ06	TA
30386056	Φ4.5×56	HAQ06	TA
30386058	Φ4.5×58	HAQ06	TA
30386060	Φ4.5×60	HAQ06	TA
30386062	Φ4.5×62	HAQ06	TA
30386064	Φ4.5×64	HAQ06	TA
30386065	Φ4.5×65	HAQ06	TA
30386070	Φ4.5×70	HAQ06	TA
30386075	Φ4.5×75	HAQ06	TA
30386080	Φ4.5×80	HAQ06	TA

SCREWS

Locking Screw, Self-tapping, 2.4



Specification	
Thread dia.	Φ2.4 mm
Shaft dia.	Φ1.9 mm
Drill bit	Φ1.8 mm

Code	Size	Model No.	Material
34270006	Φ2.4×6	SDLD03	TA
34270007	Φ2.4×7	SDLD03	TA
34270008	Φ2.4×8	SDLD03	TA
34270009	Φ2.4×9	SDLD03	TA
34270010	Φ2.4×10	SDLD03	TA
34270011	Φ2.4×11	SDLD03	TA
34270012	Φ2.4×12	SDLD03	TA
34270013	Φ2.4×13	SDLD03	TA
34270014	Φ2.4×14	SDLD03	TA
34270016	Φ2.4×16	SDLD03	TA
34270018	Φ2.4×18	SDLD03	TA
34270020	Φ2.4×20	SDLD03	TA
34270022	Φ2.4×22	SDLD03	TA
34270024	Φ2.4×24	SDLD03	TA
34270026	Φ2.4×26	SDLD03	TA
34270028	Φ2.4×28	SDLD03	TA

Code	Size	Model No.	Material
34270030	Φ2.4×30	SDLD03	TA
34270032	Φ2.4×32	SDLD03	TA
34270034	Φ2.4×34	SDLD03	TA
34270036	Φ2.4×36	SDLD03	TA
34270038	Φ2.4×38	SDLD03	TA
34270040	Φ2.4×40	SDLD03	TA
34270042	Φ2.4×42	SDLD03	TA
34270044	Φ2.4×44	SDLD03	TA
34270046	Φ2.4×46	SDLD03	TA
34270048	Φ2.4×48	SDLD03	TA
34270050	Φ2.4×50	SDLD03	TA
34270052	Φ2.4×52	SDLD03	TA
34270054	Φ2.4×54	SDLD03	TA
34270056	Φ2.4×56	SDLD03	TA
34270058	Φ2.4×58	SDLD03	TA
34270060	Φ2.4×60	SDLD03	TA

SCREWS

Locking Screw, Self-tapping, 2.7



Specification	
Thread dia.	Φ2.7 mm
Shaft dia.	Φ2.1 mm
Drill bit	Φ2.0 mm

Code	Size	Model No.	Material
34253006	Φ2.7×6	SDLD03	TA
34253007	Φ2.7×7	SDLD03	TA
34253008	Φ2.7×8	SDLD03	TA
34253009	Φ2.7×9	SDLD03	TA
34253010	Φ2.7×10	SDLD03	TA
34253011	Φ2.7×11	SDLD03	TA
34253012	Φ2.7×12	SDLD03	TA
34253013	Φ2.7×13	SDLD03	TA
34253014	Φ2.7×14	SDLD03	TA
34253016	Φ2.7×16	SDLD03	TA
34253018	Φ2.7×18	SDLD03	TA

Code	Size	Model No.	Material
34253020	Φ2.7×20	SDLD03	TA
34253022	Φ2.7×22	SDLD03	TA
34253024	Φ2.7×24	SDLD03	TA
34253026	Φ2.7×26	SDLD03	TA
34253028	Φ2.7×28	SDLD03	TA
34253030	Φ2.7×30	SDLD03	TA
34253032	Φ2.7×32	SDLD03	TA
34253034	Φ2.7×34	SDLD03	TA
34253036	Φ2.7×36	SDLD03	TA
34253038	Φ2.7×38	SDLD03	TA
34253040	Φ2.7×40	SDLD03	TA

SCREWS

Locking Screw, Self-tapping, 3.5



Specification	
Thread dia.	Φ3.5 mm
Shaft dia.	Φ2.7 mm
Drill bit	Φ2.7 mm

Code	Size	Model No.	Material
34254010	Φ3.5×10	SDLD03	TA
34254012	Φ3.5×12	SDLD03	TA
34254014	Φ3.5×14	SDLD03	TA
34254015	Φ3.5×15	SDLD03	TA
34254016	Φ3.5×16	SDLD03	TA
34254018	Φ3.5×18	SDLD03	TA
34254020	Φ3.5×20	SDLD03	TA
34254022	Φ3.5×22	SDLD03	TA
34254024	Φ3.5×24	SDLD03	TA
34254025	Φ3.5×25	SDLD03	TA
34254026	Φ3.5×26	SDLD03	TA
34254028	Φ3.5×28	SDLD03	TA
34254030	Φ3.5×30	SDLD03	TA
34254032	Φ3.5×32	SDLD03	TA

Code	Size	Model No.	Material
34254034	Φ3.5×34	SDLD03	TA
34254035	Φ3.5×35	SDLD03	TA
34254036	Φ3.5×36	SDLD03	TA
34254038	Φ3.5×38	SDLD03	TA
34254040	Φ3.5×40	SDLD03	TA
34254045	Φ3.5×45	SDLD03	TA
34254050	Φ3.5×50	SDLD03	TA
34254055	Φ3.5×55	SDLD03	TA
34254060	Φ3.5×60	SDLD03	TA
34254065	Φ3.5×65	SDLD03	TA
34254070	Φ3.5×70	SDLD03	TA
34254075	Φ3.5×75	SDLD03	TA
34254080	Φ3.5×80	SDLD03	TA

SCREWS

Locking Screw, Self-tapping, 5.0



Specification	
Thread dia.	Φ5.0 mm
Shaft dia.	Φ4.2 mm
Drill bit	Φ4.2 mm

Code	Size	Model No.	Material
34255014	Φ5.0×14	SDLD03	TA
34255015	Φ5.0×15	SDLD03	TA
34255016	Φ5.0×16	SDLD03	TA
34255018	Φ5.0×18	SDLD03	TA
34255020	Φ5.0×20	SDLD03	TA
34255022	Φ5.0×22	SDLD03	TA
34255024	Φ5.0×24	SDLD03	TA
34255025	Φ5.0×25	SDLD03	TA
34255026	Φ5.0×26	SDLD03	TA
34255028	Φ5.0×28	SDLD03	TA
34255030	Φ5.0×30	SDLD03	TA
34255032	Φ5.0×32	SDLD03	TA
34255034	Φ5.0×34	SDLD03	TA
34255035	Φ5.0×35	SDLD03	TA
34255036	Φ5.0×36	SDLD03	TA
34255038	Φ5.0×38	SDLD03	TA

Code	Size	Model No.	Material
34255040	Φ5.0×40	SDLD03	TA
34255042	Φ5.0×42	SDLD03	TA
34255044	Φ5.0×44	SDLD03	TA
34255045	Φ5.0×45	SDLD03	TA
34255046	Φ5.0×46	SDLD03	TA
34255048	Φ5.0×48	SDLD03	TA
34255050	Φ5.0×50	SDLD03	TA
34255055	Φ5.0×55	SDLD03	TA
34255060	Φ5.0×60	SDLD03	TA
34255065	Φ5.0×65	SDLD03	TA
34255070	Φ5.0×70	SDLD03	TA
34255075	Φ5.0×75	SDLD03	TA
34255080	Φ5.0×80	SDLD03	TA
34255085	Φ5.0×85	SDLD03	TA
34255090	Φ5.0×90	SDLD03	TA

SCREWS

Variable Angle Locking Screws, Self-tapping, 2.4



Specification	
Thread dia.	Φ2.4 mm
Shaft dia.	Φ1.9 mm
Drill bit	Φ1.8 mm

Code	Size	Model No.	Material
31480006	Φ2.4×6	DJY-IV	TA
31480008	Φ2.4×8	DJY-IV	TA
31480010	Φ2.4×10	DJY-IV	TA
31480012	Φ2.4×12	DJY-IV	TA
31480014	Φ2.4×14	DJY-IV	TA
31480016	Φ2.4×16	DJY-IV	TA
31480018	Φ2.4×18	DJY-IV	TA
31480020	Φ2.4×20	DJY-IV	TA
31480022	Φ2.4×22	DJY-IV	TA

Code	Size	Model No.	Material
31480024	Φ2.4×24	DJY-IV	TA
31480026	Φ2.4×26	DJY-IV	TA
31480028	Φ2.4×28	DJY-IV	TA
31480030	Φ2.4×30	DJY-IV	TA
31480032	Φ2.4×32	DJY-IV	TA
31480034	Φ2.4×34	DJY-IV	TA
31480036	Φ2.4×36	DJY-IV	TA
31480038	Φ2.4×38	DJY-IV	TA
31480040	Φ2.4×40	DJY-IV	TA

Variable Angle Locking Screws, Self-tapping, 2.7



Specification	
Thread dia.	Φ2.7 mm
Shaft dia.	Φ2.1 mm
Drill bit	Φ2.0 mm

Code	Size	Model No.	Material
31477106	Φ2.7×6	DJY-IV	TA
31477108	Φ2.7×8	DJY-IV	TA
31477110	Φ2.7×10	DJY-IV	TA
31477112	Φ2.7×12	DJY-IV	TA
31477114	Φ2.7×14	DJY-IV	TA
31477116	Φ2.7×16	DJY-IV	TA
31477118	Φ2.7×18	DJY-IV	TA
31477120	Φ2.7×20	DJY-IV	TA
31477122	Φ2.7×22	DJY-IV	TA

Code	Size	Model No.	Material
31477124	Φ2.7×24	DJY-IV	TA
31477126	Φ2.7×26	DJY-IV	TA
31477128	Φ2.7×28	DJY-IV	TA
31477130	Φ2.7×30	DJY-IV	TA
31477132	Φ2.7×32	DJY-IV	TA
31477134	Φ2.7×34	DJY-IV	TA
31477136	Φ2.7×36	DJY-IV	TA
31477138	Φ2.7×38	DJY-IV	TA
31477140	Φ2.7×40	DJY-IV	TA

SCREWS

Variable Angle Locking Screws, Self-tapping, 3.5



Specification	
Thread dia.	Φ3.5 mm
Shaft dia.	Φ2.7 mm
Drill bit	Φ2.7 mm

Code	Size	Model No.	Material
31478110	Φ3.5×10	DJY-IV	TA
31478112	Φ3.5×12	DJY-IV	TA
31478114	Φ3.5×14	DJY-IV	TA
31478116	Φ3.5×16	DJY-IV	TA
31478118	Φ3.5×18	DJY-IV	TA
31478120	Φ3.5×20	DJY-IV	TA
31478122	Φ3.5×22	DJY-IV	TA
31478124	Φ3.5×24	DJY-IV	TA
31478126	Φ3.5×26	DJY-IV	TA
31478128	Φ3.5×28	DJY-IV	TA
31478130	Φ3.5×30	DJY-IV	TA
31478132	Φ3.5×32	DJY-IV	TA

Code	Size	Model No.	Material
31478134	Φ3.5×34	DJY-IV	TA
31478136	Φ3.5×36	DJY-IV	TA
31478138	Φ3.5×38	DJY-IV	TA
31478140	Φ3.5×40	DJY-IV	TA
31478146	Φ3.5×46	DJY-IV	TA
31478150	Φ3.5×50	DJY-IV	TA
31478156	Φ3.5×56	DJY-IV	TA
31478160	Φ3.5×60	DJY-IV	TA
31478166	Φ3.5×66	DJY-IV	TA
31478170	Φ3.5×70	DJY-IV	TA
31478176	Φ3.5×76	DJY-IV	TA
31478180	Φ3.5×80	DJY-IV	TA

SCREWS

Variable Angle Locking Screws, Self-tapping, 5.0



Specification	
Thread dia.	Φ5.0 mm
Shaft dia.	Φ4.2 mm
Drill bit	Φ4.2 mm

Code	Size	Model No.	Material
31479112	Φ5.0×12	DJY-IV	TA
31479115	Φ5.0×15	DJY-IV	TA
31479118	Φ5.0×18	DJY-IV	TA
31479120	Φ5.0×20	DJY-IV	TA
31479122	Φ5.0×22	DJY-IV	TA
31479125	Φ5.0×25	DJY-IV	TA
31479128	Φ5.0×28	DJY-IV	TA
31479130	Φ5.0×30	DJY-IV	TA
31479132	Φ5.0×32	DJY-IV	TA
31479135	Φ5.0×35	DJY-IV	TA
31479138	Φ5.0×38	DJY-IV	TA
31479140	Φ5.0×40	DJY-IV	TA
31479142	Φ5.0×42	DJY-IV	TA

Code	Size	Model No.	Material
31479145	Φ5.0×45	DJY-IV	TA
31479148	Φ5.0×48	DJY-IV	TA
31479150	Φ5.0×50	DJY-IV	TA
31479155	Φ5.0×55	DJY-IV	TA
31479160	Φ5.0×60	DJY-IV	TA
31479165	Φ5.0×65	DJY-IV	TA
31479170	Φ5.0×70	DJY-IV	TA
31479175	Φ5.0×75	DJY-IV	TA
31479180	Φ5.0×80	DJY-IV	TA
31479185	Φ5.0×85	DJY-IV	TA
31479190	Φ5.0×90	DJY-IV	TA
31479195	Φ5.0×95	DJY-IV	TA
31479196	Φ5.0×100	DJY-IV	TA

SCREWS

Cannulated Locking Screw, Partial Threaded, 4.0



Specification	
Thread dia.	Φ4.0 mm
Shaft dia.	Φ3.0 mm
Inner hole	Φ1.1 mm
Drill bit size	Cannulated Drill Bit, Φ3.0/3.5mm
K wire size	Φ1.0 mm
Thread length (mm)	16/32mm

Code	Size	Model No.	Material
31475232	Φ4.0×32	DJY-IV	TA
31475240	Φ4.0×40	DJY-IV	TA
31475248	Φ4.0×48	DJY-IV	TA
31475255	Φ4.0×55	DJY-IV	TA
31475258	Φ4.0×58	DJY-IV	TA
31475262	Φ4.0×62	DJY-IV	TA
31475265	Φ4.0×65	DJY-IV	TA

Code	Size	Model No.	Material
31475270	Φ4.0×70	DJY-IV	TA
31475272	Φ4.0×72	DJY-IV	TA
31475278	Φ4.0×78	DJY-IV	TA
31475280	Φ4.0×80	DJY-IV	TA
31475285	Φ4.0×85	DJY-IV	TA
31475288	Φ4.0×88	DJY-IV	TA

Cannulated Locking Screw, Partial Threaded, 5.5



Specification	
Thread dia.	Φ5.5 mm
Shaft dia.	Φ4.5 mm
Inner hole	Φ1.7 mm
Drill bit size	Cannulated Drill Bit, Φ4.5/5.0mm
K wire size	Φ1.6 mm
Thread length (mm)	16/22/32mm

Code	Size	Model No.	Material
31476250	Φ5.5×50	DJY-IV	TA
31476255	Φ5.5×55	DJY-IV	TA
31476262	Φ5.5×62	DJY-IV	TA
31476265	Φ5.5×65	DJY-IV	TA
31476270	Φ5.5×70	DJY-IV	TA
31476272	Φ5.5×72	DJY-IV	TA

Code	Size	Model No.	Material
31476278	Φ5.5×78	DJY-IV	TA
31476280	Φ5.5×80	DJY-IV	TA
31476285	Φ5.5×85	DJY-IV	TA
31476292	Φ5.5×92	DJY-IV	TA
31476293	Φ5.5×100	DJY-IV	TA

SCREWS

Cannulated Locking Screw, Partial Threaded, 7.0



Specification	
Thread dia.	Φ7.0 mm
Shaft dia.	Φ4.5 mm
Inner hole	Φ2.6 mm
Drill bit size	Cannulated Drill Bit, Φ 4.5/6.5mm
K wire size	Φ2.5 mm
Thread length (mm)	32mm

Code	Size	Model No.	Material
31474070	Φ7.0×70	DJY-IV	TA
31474075	Φ7.0×75	DJY-IV	TA
31474078	Φ7.0×78	DJY-IV	TA
31474082	Φ7.0×82	DJY-IV	TA
31474085	Φ7.0×85	DJY-IV	TA
31474090	Φ7.0×90	DJY-IV	TA
31474095	Φ7.0×95	DJY-IV	TA
31474100	Φ7.0×100	DJY-IV	TA


Code	Size	Model No.	Material
31474105	Φ7.0×105	DJY-IV	TA
31474110	Φ7.0×110	DJY-IV	TA
31474115	Φ7.0×115	DJY-IV	TA
31474120	Φ7.0×120	DJY-IV	TA
31474125	Φ7.0×125	DJY-IV	TA
31474130	Φ7.0×130	DJY-IV	TA
31474135	Φ7.0×135	DJY-IV	TA
31474140	Φ7.0×140	DJY-IV	TA

INSTRUMENTS

INSTRUMENTS

VariLoc 2.4/2.7mm Plate System

Code	Description	QTY.
295000	VariLoc 2.4/2.7mm Plate System	1
295070	VariLoc 2.4/2.7mm Plate System-empty	1
295050	VariLoc 2.4/2.7mm Screw Case	1
286120	Drill Guide 2.0	1
286130	Drill Bit, Φ 2.0mm	1
292110	Drill Bit, Φ 2.0mm	1
286220	Drill Guide 1.8	2
286230	Drill Bit, Φ 1.8mm	2
295110	Drill Bit, Φ 1.8mm	2
288150	Torque-limiting Screwdriver,T8	2
288160	Locking Screw Holding Sleeve	1
288170	Cortex Screw Holding Sleeve	1
286200	Depth Gauge,60mm	1
288130	Torx Screwdriver, with Quick Coupling,T8	1
286260	Drill Bit, Φ 2.7mm	1
286270	Drill Bit, Φ2.4mm	1
286240	Double Drill Guide 2.7/2.0	1
286250	Double Drill Guide 2.4/1.8	1
292140	Drill Guide,for 2.7 VA Screw	2
295112	Drill Guide,for 2.4 VA Screw	2
10735150	K-Wire, ø1.8x150mm	3
19472515	K-Wire with thread, ø1.8x150mm	1
288140	Screwdriver for Screw Removal, T8	1
01214	Bending Iron	1
01123	Reduction Forceps with Points	1
01124	Reduction Forceps with Points, toothed	1
01232	Periosteal Elevator, width 3mm	1
01215	Retractor, width 6mm	1
291290	T-handle with Quick Coupling	1



INSTRUMENTS

INSTRUMENTS



295050
VariLoc 2.4/2.7mm Screw Case



286120
Drill Guide 2.0



286130
Drill Bit, Φ 2.0mm



292110
Drill Bit, Φ 2.0mm



286220
Drill Guide 1.8



286230
Drill Bit, Φ 1.8mm



295110
Drill Bit, Φ 1.8mm



288150
Torque-limiting Screwdriver,T8



288160
Locking Screw Holding Sleeve



288170
Cortex Screw Holding Sleeve



286200
Depth Gauge,60mm



288130
Torx Screwdriver, with Quick Coupling,T8

INSTRUMENTS

INSTRUMENTS



286260
Drill Bit, Φ 2.7mm



286270
Drill Bit, Φ2.4mm



286240
Double Drill Guide 2.7/2.0



286250
Double Drill Guide 2.4/1.8



292140
Drill Guide,for 2.7 VA Screw



295112
Drill Guide,for 2.4 VA Screw



10735150
K-Wire, ø1.8x150mm



19472515
K-Wire with thread, ø1.8x150mm



288140
Screwdriver for Screw Removal, T8



01214
Bending Iron



01123
Reduction Forceps with Points



01124
Reduction Forceps with Points, toothed

INSTRUMENTS

INSTRUMENTS




INSTRUMENTS

INSTRUMENTS

VariLoc 2.7/3.5mm Plate System

Code	Description	QTY.
292000	VariLoc 2.7/3.5mm Plate System	1
292010	VariLoc 2.7/3.5mm Plate System-empty	1
292050	VariLoc 2.7/3.5mm Screw Case	1
286130	Drill Bit, Φ 2.0mm	2
292110	Drill Bit, Φ 2.0mm	2
286260	Drill Bit, Φ 2.7mm	2
286120	Drill Sleeve for Drill Bit Φ 2.0mm	2
286240	Double Drill Guide 2.0/2.7	1
286280	Tap, ø2.7mm	1
288150	Torque-limiting Screwdriver,T8	1
288130	Torx Screwdriver, with Quick Coupling, T8	1
288140	Screwdriver for Screw Removal, T8	1
288160	Locking Screw Holding Sleeve	1
288170	Cortex Screw Holding Sleeve	1
286200	Depth Gauge	1
292580	Drill Bit with Stop, Φ 2.7mm	2
292220	Drill Bit with Stop, Φ 2.7mm	2
292240	Drill Sleeve for Drill Bit Φ 2.7mm	3
010010	Drill Bit, Φ 2.5mm	1
292230	Drill Bit, Φ 2.5mm	1
010030	Drill Bit, Φ 3.5mm	1
201120	Double Drill Guide 2.5/3.5	1
020010	Tap, ø3.5mm	1
292280	Torque-limiting Screwdriver, T15	1
292290	Torx Screwdriver, with Quick Coupling,T15	1
292300	Screwdriver for Screw Removal, T15	1
292310	Locking Screw Holding Sleeve, Φ 3.5mm	1
292320	Cortex Screw Holding Sleeve, Φ 3.5mm	1
020020	Tap, ø4.0mm	1



INSTRUMENTS

INSTRUMENTS

VariLoc 2.7/3.5mm Plate System

Code	Description	QTY.
292450	Hex Screwdriver	1
201200	Depth Gauge	1
10736220	K-Wire with thread, Φ 1.2mm	3
10703210	K-Wire, Φ 1.2mm	3
221180	Allen Key	1
291290	T-handle with Quick Coupling	1
201170	Bone Holding Forceps, self-centering	2
201180	Reduction Forceps, serrated jaws	1
201190	Reduction Forceps with Points	1
01214	Bending Iron, for 2.7mm LCP	2
201160	Bending Iron, for 3.5mm LCP	2
206270	Bending Iron, for Recon LCP	2
10702200	K-Wire, Φ 1.0mm	3
291200	K-Wire with thread, Φ 1.0mm	3
292360	Guide Sleeve for K-Wire, Φ 1.0mm	1
292370	Depth Gauge for K-Wire, Φ 1.0mm	1
292390	Drill Guide ø3.5	2
292610	Cannulated Drill Bit, Φ3.0/3.5mm	2
292680	Cannulated Drill Bit, Φ3.0/3.5mm	2
292410	Cannulated Torx Screwdriver, with Quick Coupling,T15	1
292140	Drill Guide,for 2.7 VA Screw	1
292250	Drill Guide,for 3.5 VA Screw	1





292050
VariLoc 2.7/3.5mm Screw Case



286130
Drill Bit, Φ 2.0mm



292110
Drill Bit, Φ 2.0mm

INSTRUMENTS

INSTRUMENTS



286260
Drill Bit, Φ 2.7mm



286120
Drill Sleeve for Drill Bit Φ 2.0mm



286240
Double Drill Guide 2.0/2.7



286280
Tap, ø2.7mm



288150
Torque-limiting Screwdriver,T8



288130
Torx Screwdriver, with Quick Coupling, T8



288140
Screwdriver for Screw Removal, T8



288160
Locking Screw Holding Sleeve



288170
Cortex Screw Holding Sleeve



286200
Depth Gauge



292580
Drill Bit with Stop, Φ 2.7mm



292220
Drill Bit with Stop, Φ 2.7mm

INSTRUMENTS

INSTRUMENTS

 292240 Drill Sleeve for Drill Bit Φ 2.7mm	 010010 Drill Bit, Φ 2.5mm	 292230 Drill Bit, Φ 2.5mm
 010030 Drill Bit, Φ 3.5mm	 201120 Double Drill Guide 2.5/3.5	 020010 Tap, $\varnothing$3.5mm
 292280 Torque-limiting Screwdriver, T15	 292290 Torx Screwdriver, with Quick Coupling,T15	 292300 Screwdriver for Screw Removal, T15
 292310 Locking Screw Holding Sleeve, Φ 3.5mm	 292320 Cortex Screw Holding Sleeve, Φ 3.5mm	 020020 Tap, $\varnothing$4.0mm


INSTRUMENTS

INSTRUMENTS

 292450 Hex Screwdriver	 201200 Depth Gauge	 10736220 K-Wire with thread, Φ 1.2mm
 10703210 K-Wire, Φ 1.2mm	 221180 Allen Key	 291290 T-handle with Quick Coupling
 201170 Bone Holding Forceps, self-centering	 201180 Reduction Forceps, serrated jaws	 201190 Reduction Forceps with Points
 01214 Bending Iron, for 2.7mm LCP	 201160 Bending Iron, for 3.5mm LCP	 206270 Bending Iron, for Recon LCP

INSTRUMENTS

INSTRUMENTS



291200
K-Wire with thread, Ø 1.0mm



292360
Guide Sleeve for K-Wire, Ø 1.0mm



292370
Depth Gauge for K-Wire, Ø 1.0mm



292390
Drill Guide ø3.5



292610
Cannulated Drill Bit, Ø3.0/3.5mm



292680
Cannulated Drill Bit, Ø3.0/3.5mm




292410
Cannulated Torx Screwdriver, with Quick Coupling,T15



292140
Drill Guide,for 2.7 VA Screw



292250
Drill Guide,for 3.5 VA Screw



10702200
K-Wire, 1.0mm

INSTRUMENTS

INSTRUMENTS

VariLoc 5.0mm Plate System

Code	Description	QTY.
291000	VariLoc 5.0mm Plate System	
291010	VariLoc 5.0mm Plate System-empty	1
291050	VariLoc 5.0mm Screw Case	1
291290	T-handle with Quick Coupling	1
221180	Allen Key	1
010020	Drill Bit , Ø 3.2mm	1
291100	Drill Bit , Ø 3.2mm	1
010040	Drill Bit , Ø 4.5mm	1
020030	Tap,Ø 4.5mm	1
206120	Double Drill Guide 4.5/3.2	1
291600	Drill Bit with Stop, Ø 4.2mm	2
291130	Drill Bit with Stop, Ø 4.2mm	2
291140	Drill Guide 4.2	3
291150	Torque-limiting Screwdriver, T25	1
291160	Torx Screwdriver, with Quick Coupling,T25	1
291170	Screwdriver for Screw Removal, T25	1
291270	Locking Screw Holding Sleeve, Ø 5.0mm	1
291180	Cortex Screw Holding Sleeve, Ø 4.5mm	1
020040	Tap,Ø 6.5mm	1
206110	Double Drill Guide 6.5/3.2	1
291320	Hex Screwdriver	1
206170	Depth Gauge	1
10737230	K-Wire with thread, Ø 2.0mm	3
10705230	K-Wire, Ø 2.0mm	3
206290	Bone Holding Forceps, self-centering	2
206150	Reduction Forceps, serrated jaws, large	1
206160	Reduction Forceps with Points, large	1
206320	Bending Irons(left)	1
206340	Bending Irons(right)	1
291190	K-Wire with thread, Ø 1.6mm	3
19472425	K-Wire, Ø 1.6mm	3
291210	Guide Sleeve for K-Wire, Ø 1.6mm	1
291220	Depth Gauge for K-Wire, Ø 1.6mm	1
291240	Drill Guide 5.0	2
291250	Cannulated Tap,Ø 5.5mm	1
291620	Cannulated Torx Screwdriver, with Quick Coupling,T25	1
291630	Cannulated Drill Bit, Ø4.5/5.0mm	2
291570	Cannulated Drill Bit, Ø4.5/5.0mm	2
291710	Drill Bit with Stop, Ø 4.2mm for 5.0 VA Screw	2
291720	Drill Guide ,for 5.0 VA Screw	1



INSTRUMENTS

INSTRUMENTS



291050
VariLoc 5.0mm Screw Case



291290
T-handle with Quick Coupling



221180
Allen Key



010020
Drill Bit , Ø 3.2mm



291100
Drill Bit , Ø 3.2mm



010040
Drill Bit , Ø 4.5mm



020030
Tap, Ø 4.5mm



206120
Double Drill Guide 4.5/3.2



291600
Drill Bit with Stop, Ø 4.2mm



291130
Drill Bit with Stop, Ø 4.2mm



291140
Drill Guide 4.2



291150
Torque-limiting Screwdriver, T25

INSTRUMENTS

INSTRUMENTS



291160
Torx Screwdriver,
with Quick Coupling,T25



291170
Screwdriver for Screw Removal,
T25



291270
Locking Screw Holding Sleeve,
Ø 5.0mm



291180
Cortex Screw Holding Sleeve,
Ø 4.5mm



020040
Tap, Ø 6.5mm



206110
Double Drill Guide 6.5/3.2



291320
Hex Screwdriver



206170
Depth Gauge



10737230
K-Wire with thread, Ø 2.0mm



10705230
K-Wire, Ø 2.0mm



206290
Bone Holding Forceps,
self-centering



206150
Reduction Forceps,
serrated jaws, large

INSTRUMENTS

INSTRUMENTS



206160
Reduction Forceps with Points,
large



206320
Bending Irons(left)



206340
Bending Irons(right)



291190
K-Wire with thread, Ø 1.6mm



19472425
K-Wire, Ø 1.6mm



291210
Guide Sleeve for K-Wire, Ø 1.6mm



291220
Depth Gauge for K-Wire, Ø 1.6mm



291240
Drill Guide 5.0



291250
Cannulated Tap, Ø 5.5mm



291620
Cannulated Torx Screwdriver,
with Quick Coupling,T25



291630
Cannulated Drill Bit, Ø4.5/5.0mm



291570
Cannulated Drill Bit, Ø4.5/5.0mm

INSTRUMENTS

INSTRUMENTS



291710
Drill Bit with Stop,
4.2mm for 5.0 VA Screw




291720
Drill Guide ,for 5.0 VA Screw

INSTRUMENTS

INSTRUMENTS


VariLoc 5.0/7.0mm Proximal Femoral Plate System

Code	Description	QTY.
293000	VariLoc 5.0/7.0mm Proximal Femoral Plate System	1
293010	VariLoc 5.0/7.0mm Proximal Femoral Plate System-empty case	1
293050	5.0/7.0mm Screw Case	1
10737230	K-Wire with thread, Φ 2.0mm	3
10705230	K-Wire, Φ 2.0mm	3
291740	K-Wire, Φ 2.5mm	3
291750	K-Wire with thread, Φ 2.5mm	3
291410	Guide Sleeve for K-Wire, Φ 2.5mm	1
291400	Depth Gauge for K-Wire, Φ 2.5mm	1
291430	Drill Sleeve,Φ 6.5mm	3
291260	Cannulated Drill Bit, Φ 4.5/6.5mm	2
291440	Cannulated Drill Bit, Φ 4.5/6.5mm	2
291450	Cannulated Tap,Φ 7.0mm	1
291670	Drill Sleeve for Cannulated Drill Bit	2
291470	Torque-limiting Screwdriver,Φ 7.0mm	1
291480	Torx Screwdriver, Φ 7.0mm with Quick Coupling	1
291490	Screwdriver for Screw Removal	1
291500	Screw Holding Sleeve, Φ 7.0mm	1
291540	T-handle with Quick Coupling	1
291600	Drill Bit with Stop, Φ 4.2mm	2
291130	Drill Bit with Stop, Φ 4.2mm	2
291140	Drill Guide 4.2	2
291150	Torque-limiting Screwdriver, T25	1
291160	Torx Screwdriver, with Quick Coupling,T25	1
291170	Screwdriver for Screw Removal, T25	1
291270	Locking Screw Holding Sleeve, Φ 5.0mm	1
221180	Allen Key	1
291290	T-handle with Quick Coupling	1



INSTRUMENTS

INSTRUMENTS

 293050 5.0/7.0mm Screw Case	 10737230 K-Wire with thread, Φ 2.0mm	 10705230 K-Wire, Φ 2.0mm
 291740 K-Wire, Φ 2.5mm	 291750 K-Wire with thread, Φ 2.5mm	 291410 Guide Sleeve for K-Wire, Φ 2.5mm
 291400 Depth Gauge for K-Wire, Φ 2.5mm	 291430 Drill Sleeve, Φ 6.5mm	 291260 Cannulated Drill Bit, Φ 4.5/6.5mm
 291440 Cannulated Drill Bit, Φ 4.5/6.5mm	 291450 Cannulated Tap, Φ 7.0mm	 291670 Drill Sleeve for Cannulated Drill Bit

INSTRUMENTS

INSTRUMENTS



291470
Torque-limiting Screwdriver,
Φ 7.0mm



291480
Torx Screwdriver, Φ 7.0mm with
Quick Coupling



291490
Screwdriver for Screw Removal




291500
Screw Holding Sleeve, Φ 7.0mm



291540
T-handle with Quick Coupling



291600
Drill Bit with Stop, Φ 4.2mm



291130
Drill Bit with Stop, Φ 4.2mm



291140
Drill Guide 4.2



291150
Torque-limiting Screwdriver, T25



291160
Torx Screwdriver, with Quick
Coupling,T25



291170
Screwdriver for Screw Removal,
T25



291270
Locking Screw Holding Sleeve,
Φ 5.0mm

INSTRUMENTS

INSTRUMENTS



221180
Allen Key




291290
T-handle with Quick Coupling

INSTRUMENTS

INSTRUMENTS

VariLoc Locking Compression Plates Removal Kit

Code	Description	QTY.
294000	VariLoc LCP Removal Kit	1
294010	VariLoc LCP Removal Kit-empty case	1
288120	Easyout for Screw ø2.7	1
292330	Easyout for Screw ø3.5	1
291360	Easyout for Screw ø5.0	1
291530	Easyout for Screw ø7.0	1
221190	Drill Bit for Screw Head Removal,ø2.5	1
226190	Drill Bit for Screw Head Removal,ø3.5	1
292670	Trephine, ø3.5	1
291350	Trephine, ø5.0	1
291290	T-handle with Quick Couping	1
288140	Screwdriver for Screw Removal, T8	1
292300	Screwdriver for Screw Removal, T15	1
291170	Screwdriver for Screw Removal, T25	1
291490	Screwdriver for Screw Removal,ø7.0	1
292450	Hex Screwdriver SW2.5	1
291320	Hex Screwdriver SW3.5	1





291290
T-handle with Quick Couping



288140
Screwdriver for Screw Removal, T8



291170
Screwdriver for Screw Removal, T25

INSTRUMENTS

INSTRUMENTS



288120
Easyout for Screw ø2.7



292330
Easyout for Screw ø3.5



291360
Easyout for Screw ø5.0



291530
Easyout for Screw ø7.0



221190
Drill Bit for Screw Head Removal,ø2.5



226190
Drill Bit for Screw Head Removal,ø3.5



292670
Trephine, ø3.5



291350
Trephine, ø5.0



291320
Hex Screwdriver SW3.5



292300
Screwdriver for Screw Removal, T15



291490
Screwdriver for Screw Removal, ø7.0



292450
Hex Screwdriver SW2.5

INSTRUMENTS

SCREW CASE

295050 VariLoc 2.4/2.7mm Screw Case



Code	Description	Size	QTY.
34270010	Locking Screw, self-tapping	Ø2.4×10	12
34270012	Locking Screw, self-tapping	Ø2.4×12	12
34270014	Locking Screw, self-tapping	Ø2.4×14	12
34270016	Locking Screw, self-tapping	Ø2.4×16	12
34270018	Locking Screw, self-tapping	Ø2.4×18	12
34270020	Locking Screw, self-tapping	Ø2.4×20	12
34270022	Locking Screw, self-tapping	Ø2.4×22	12
34270024	Locking Screw, self-tapping	Ø2.4×24	12
34270026	Locking Screw, self-tapping	Ø2.4×26	12
34270028	Locking Screw, self-tapping	Ø2.4×28	12
34270030	Locking Screw, self-tapping	Ø2.4×30	12
30374010	Cortex Screw, fully threaded,self-tapping	Ø2.4×10	6
30374012	Cortex Screw, fully threaded,self-tapping	Ø2.4×12	6
30374014	Cortex Screw, fully threaded,self-tapping	Ø2.4×14	6
30374016	Cortex Screw, fully threaded,self-tapping	Ø2.4×16	6
30374018	Cortex Screw, fully threaded,self-tapping	Ø2.4×18	6
30374020	Cortex Screw, fully threaded,self-tapping	Ø2.4×20	6
30374022	Cortex Screw, fully threaded,self-tapping	Ø2.4×22	6
30374024	Cortex Screw, fully threaded,self-tapping	Ø2.4×24	6
30374026	Cortex Screw, fully threaded,self-tapping	Ø2.4×26	6
30374028	Cortex Screw, fully threaded,self-tapping	Ø2.4×28	6
30374030	Cortex Screw, fully threaded,self-tapping	Ø2.4×30	6
31480010	Variable Angle Locking Screws, self-tapping	Ø2.4×10	4
31480012	Variable Angle Locking Screws, self-tapping	Ø2.4×12	4
31480014	Variable Angle Locking Screws, self-tapping	Ø2.4×14	4
31480016	Variable Angle Locking Screws, self-tapping	Ø2.4×16	4
31480018	Variable Angle Locking Screws, self-tapping	Ø2.4×18	4
31480020	Variable Angle Locking Screws, self-tapping	Ø2.4×20	4
31480022	Variable Angle Locking Screws, self-tapping	Ø2.4×22	4
31480024	Variable Angle Locking Screws, self-tapping	Ø2.4×24	4
31480026	Variable Angle Locking Screws, self-tapping	Ø2.4×26	4
31480028	Variable Angle Locking Screws, self-tapping	Ø2.4×28	4
31480030	Variable Angle Locking Screws, self-tapping	Ø2.4×30	4
34253010	Locking Screw, self-tapping	Ø2.7×10	12
34253012	Locking Screw, self-tapping	Ø2.7×12	12
34253014	Locking Screw, self-tapping	Ø2.7×14	12
34253016	Locking Screw, self-tapping	Ø2.7×16	12
34253018	Locking Screw, self-tapping	Ø2.7×18	12
34253020	Locking Screw, self-tapping	Ø2.7×20	12
34253022	Locking Screw, self-tapping	Ø2.7×22	12
34253024	Locking Screw, self-tapping	Ø2.7×24	12

INSTRUMENTS

SCREW CASE

Code	Description	Size	QTY.
34253026	Locking Screw, self-tapping	Ø2.7×26	12
34253028	Locking Screw, self-tapping	Ø2.7×28	12
34253030	Locking Screw, self-tapping	Ø2.7×30	12
34253032	Locking Screw, self-tapping	Ø2.7×32	12
34253034	Locking Screw, self-tapping	Ø2.7×34	12
34253036	Locking Screw, self-tapping	Ø2.7×36	12
34253038	Locking Screw, self-tapping	Ø2.7×38	12
34253040	Locking Screw, self-tapping	Ø2.7×40	12
30375010	Cortex Screw, fully threaded,self-tapping	Ø2.7×10	6
30375012	Cortex Screw, fully threaded,self-tapping	Ø2.7×12	6
30375014	Cortex Screw, fully threaded,self-tapping	Ø2.7×14	6
30375016	Cortex Screw, fully threaded,self-tapping	Ø2.7×16	6
30375018	Cortex Screw, fully threaded,self-tapping	Ø2.7×18	6
30375020	Cortex Screw, fully threaded,self-tapping	Ø2.7×20	6
30375022	Cortex Screw, fully threaded,self-tapping	Ø2.7×22	6
30375024	Cortex Screw, fully threaded,self-tapping	Ø2.7×24	6
30375026	Cortex Screw, fully threaded,self-tapping	Ø2.7×26	6
30375028	Cortex Screw, fully threaded,self-tapping	Ø2.7×28	6
30375030	Cortex Screw, fully threaded,self-tapping	Ø2.7×30	6
30375032	Cortex Screw, fully threaded,self-tapping	Ø2.7×32	6
30375034	Cortex Screw, fully threaded,self-tapping	Ø2.7×34	6
30375036	Cortex Screw, fully threaded,self-tapping	Ø2.7×36	6
30375038	Cortex Screw, fully threaded,self-tapping	Ø2.7×38	6
30375040	Cortex Screw, fully threaded,self-tapping	Ø2.7×40	6
31477110	Variable Angle Locking Screws, self-tapping	Ø2.7×10	4
31477112	Variable Angle Locking Screws, self-tapping	Ø2.7×12	4
31477114	Variable Angle Locking Screws, self-tapping	Ø2.7×14	4
31477116	Variable Angle Locking Screws, self-tapping	Ø2.7×16	4
31477118	Variable Angle Locking Screws, self-tapping	Ø2.7×18	4
31477120	Variable Angle Locking Screws, self-tapping	Ø2.7×20	4
31477122	Variable Angle Locking Screws, self-tapping	Ø2.7×22	4
31477124	Variable Angle Locking Screws, self-tapping	Ø2.7×24	4
31477126	Variable Angle Locking Screws, self-tapping	Ø2.7×26	4
31477128	Variable Angle Locking Screws, self-tapping	Ø2.7×28	4
31477130	Variable Angle Locking Screws, self-tapping	Ø2.7×30	4
31477132	Variable Angle Locking Screws, self-tapping	Ø2.7×32	4
31477134	Variable Angle Locking Screws, self-tapping	Ø2.7×34	4
31477136	Variable Angle Locking Screws, self-tapping	Ø2.7×36	4
31477138	Variable Angle Locking Screws, self-tapping	Ø2.7×38	4
31477140	Variable Angle Locking Screws, self-tapping	Ø2.7×40	4

INSTRUMENTS

SCREW CASE

292050 VariLoc 2.7/3.5mm Screw Case



Code	Description	Size	QTY.
34253010	Locking Screw, self-tapping	Ø2.7×10	6
34253012	Locking Screw, self-tapping	Ø2.7×12	6
34253014	Locking Screw, self-tapping	Ø2.7×14	6
34253016	Locking Screw, self-tapping	Ø2.7×16	6
34253018	Locking Screw, self-tapping	Ø2.7×18	6
34253020	Locking Screw, self-tapping	Ø2.7×20	6
34253022	Locking Screw, self-tapping	Ø2.7×22	6
34253024	Locking Screw, self-tapping	Ø2.7×24	6
34253026	Locking Screw, self-tapping	Ø2.7×26	6
34253028	Locking Screw, self-tapping	Ø2.7×28	6
34253030	Locking Screw, self-tapping	Ø2.7×30	6
34253032	Locking Screw, self-tapping	Ø2.7×32	6
34253034	Locking Screw, self-tapping	Ø2.7×34	6
34253036	Locking Screw, self-tapping	Ø2.7×36	6
34253038	Locking Screw, self-tapping	Ø2.7×38	6
34253040	Locking Screw, self-tapping	Ø2.7×40	6
30375010	Cortex Screw, fully threaded,self-tapping	Ø2.7×10	3
30375012	Cortex Screw, fully threaded,self-tapping	Ø2.7×12	3
30375014	Cortex Screw, fully threaded,self-tapping	Ø2.7×14	3
30375016	Cortex Screw, fully threaded,self-tapping	Ø2.7×16	3
30375018	Cortex Screw, fully threaded,self-tapping	Ø2.7×18	3
30375020	Cortex Screw, fully threaded,self-tapping	Ø2.7×20	3
30375022	Cortex Screw, fully threaded,self-tapping	Ø2.7×22	3
30375024	Cortex Screw, fully threaded,self-tapping	Ø2.7×24	3
30375026	Cortex Screw, fully threaded,self-tapping	Ø2.7×26	3
30375028	Cortex Screw, fully threaded,self-tapping	Ø2.7×28	3
30375030	Cortex Screw, fully threaded,self-tapping	Ø2.7×30	3
30375032	Cortex Screw, fully threaded,self-tapping	Ø2.7×32	3
30375034	Cortex Screw, fully threaded,self-tapping	Ø2.7×34	3
30375036	Cortex Screw, fully threaded,self-tapping	Ø2.7×36	3
30375038	Cortex Screw, fully threaded,self-tapping	Ø2.7×38	3
30375040	Cortex Screw, fully threaded,self-tapping	Ø2.7×40	3
34254010	Locking Screw, self-tapping	Ø3.5×10	6
34254012	Locking Screw, self-tapping	Ø3.5×12	6
34254014	Locking Screw, self-tapping	Ø3.5×14	6
34254016	Locking Screw, self-tapping	Ø3.5×16	6
34254018	Locking Screw, self-tapping	Ø3.5×18	6
34254020	Locking Screw, self-tapping	Ø3.5×20	6
34254022	Locking Screw, self-tapping	Ø3.5×22	6
34254024	Locking Screw, self-tapping	Ø3.5×24	6
34254026	Locking Screw, self-tapping	Ø3.5×26	6

INSTRUMENTS

SCREW CASE

Code	Description	Size	QTY.
34254028	Locking Screw, self-tapping	Ø3.5×28	6
34254030	Locking Screw, self-tapping	Ø3.5×30	6
34254032	Locking Screw, self-tapping	Ø3.5×32	6
34254034	Locking Screw, self-tapping	Ø3.5×34	6
34254036	Locking Screw, self-tapping	Ø3.5×36	6
34254038	Locking Screw, self-tapping	Ø3.5×38	6
34254040	Locking Screw, self-tapping	Ø3.5×40	6
34254045	Locking Screw, self-tapping	Ø3.5×45	6
34254050	Locking Screw, self-tapping	Ø3.5×50	6
34254055	Locking Screw, self-tapping	Ø3.5×55	6
34254060	Locking Screw, self-tapping	Ø3.5×60	6
34254065	Locking Screw, self-tapping	Ø3.5×65	6
34254070	Locking Screw, self-tapping	Ø3.5×70	6
30387010	Cortex Screw, fully threaded,self-tapping	Ø3.5x10	3
30387012	Cortex Screw, fully threaded,self-tapping	Ø3.5x12	3
30387014	Cortex Screw, fully threaded,self-tapping	Ø3.5x14	3
30387016	Cortex Screw, fully threaded,self-tapping	Ø3.5x16	3
30387018	Cortex Screw, fully threaded,self-tapping	Ø3.5x18	3
30387020	Cortex Screw, fully threaded,self-tapping	Ø3.5x20	3
30387022	Cortex Screw, fully threaded,self-tapping	Ø3.5x22	3
30387024	Cortex Screw, fully threaded,self-tapping	Ø3.5x24	3
30387026	Cortex Screw, fully threaded,self-tapping	Ø3.5x26	3
30387028	Cortex Screw, fully threaded,self-tapping	Ø3.5x28	3
30387030	Cortex Screw, fully threaded,self-tapping	Ø3.5x30	3
30387032	Cortex Screw, fully threaded,self-tapping	Ø3.5x32	3
30387034	Cortex Screw, fully threaded,self-tapping	Ø3.5x34	3
30387036	Cortex Screw, fully threaded,self-tapping	Ø3.5x36	3
30387038	Cortex Screw, fully threaded,self-tapping	Ø3.5x38	3
30387040	Cortex Screw, fully threaded,self-tapping	Ø3.5x40	3
30387045	Cortex Screw, fully threaded,self-tapping	Ø3.5x45	3
30387050	Cortex Screw, fully threaded,self-tapping	Ø3.5x50	3
30387055	Cortex Screw, fully threaded,self-tapping	Ø3.5x55	3
30387060	Cortex Screw, fully threaded,self-tapping	Ø3.5x60	3
30403110	Cancellous Screws, fully threaded	4.0×10mm	2
30403112	Cancellous Screws, fully threaded	4.0×12mm	2
30403114	Cancellous Screws, fully threaded	4.0×14mm	2
30403116	Cancellous Screws, fully threaded	4.0×16mm	2
30403118	Cancellous Screws, fully threaded	4.0×18mm	2
30403120	Cancellous Screws, fully threaded	4.0×20mm	2
30403122	Cancellous Screws, fully threaded	4.0×22mm	2
30403124	Cancellous Screws, fully threaded	4.0×24mm	2

INSTRUMENTS

SCREW CASE

Code	Description	Size	QTY.
30403126	Cancellous Screws, fully threaded	4.0×26mm	2
30403128	Cancellous Screws, fully threaded	4.0×28mm	2
30403130	Cancellous Screws, fully threaded	4.0×30mm	2
30403132	Cancellous Screws, fully threaded	4.0×32mm	2
30403134	Cancellous Screws, fully threaded	4.0×34mm	2
30403135	Cancellous Screws, fully threaded	4.0×36mm	2
30403136	Cancellous Screws, fully threaded	4.0×38mm	2
30403138	Cancellous Screws, fully threaded	4.0×40mm	2
30403140	Cancellous Screws, fully threaded	4.0×42mm	2
30403142	Cancellous Screws, fully threaded	4.0×44mm	2
30403144	Cancellous Screws, fully threaded	4.0×45mm	2
30403145	Cancellous Screws, fully threaded	4.0×46mm	2
30403146	Cancellous Screws, fully threaded	4.0×48mm	2
30403148	Cancellous Screws, fully threaded	4.0×50mm	2
30403150	Cancellous Screws, fully threaded	4.0×52mm	2
30403152	Cancellous Screws, fully threaded	4.0×55mm	2
30403155	Cancellous Screws, fully threaded	4.0×56mm	2
30403160	Cancellous Screws, fully threaded	4.0×60mm	2
30403165	Cancellous Screws, fully threaded	4.0×65mm	2
30403170	Cancellous Screws, fully threaded	4.0×70mm	2
30403210	Cancellous Screws, partially threaded	4.0×10mm	2
30403212	Cancellous Screws, partially threaded	4.0×12mm	2
30403214	Cancellous Screws, partially threaded	4.0×14mm	2
30403216	Cancellous Screws, partially threaded	4.0×16mm	2
30403218	Cancellous Screws, partially threaded	4.0×18mm	2
30403220	Cancellous Screws, partially threaded	4.0×20mm	2
30403222	Cancellous Screws, partially threaded	4.0×22mm	2
30403224	Cancellous Screws, partially threaded	4.0×24mm	2
30403226	Cancellous Screws, partially threaded	4.0×26mm	2
30403228	Cancellous Screws, partially threaded	4.0×28mm	2
30403230	Cancellous Screws, partially threaded	4.0×30mm	2
30403232	Cancellous Screws, partially threaded	4.0×32mm	2
30403234	Cancellous Screws, partially threaded	4.0×34mm	2
30403236	Cancellous Screws, partially threaded	4.0×36mm	2
30403238	Cancellous Screws, partially threaded	4.0×38mm	2
30403240	Cancellous Screws, partially threaded	4.0×40mm	2
30403242	Cancellous Screws, partially threaded	4.0×42mm	2
30403244	Cancellous Screws, partially threaded	4.0×44mm	2
30403246	Cancellous Screws, partially threaded	4.0×46mm	2
30403248	Cancellous Screws, partially threaded	4.0×48mm	2
30403250	Cancellous Screws, partially threaded	4.0×50mm	2

INSTRUMENTS

SCREW CASE

Code	Description	Size	QTY.
30403252	Cancellous Screws, partially threaded	4.0×52mm	2
30403254	Cancellous Screws, partially threaded	4.0×54mm	2
30403256	Cancellous Screws, partially threaded	4.0×56mm	2
30403260	Cancellous Screws, partially threaded	4.0×60mm	2
30403265	Cancellous Screws, partially threaded	4.0×65mm	2
30403270	Cancellous Screws, partially threaded	4.0×70mm	2
31477110	Variable Angle Locking Screws, self-tapping	Φ2.7×10	2
31477112	Variable Angle Locking Screws, self-tapping	Φ2.7×12	2
31477114	Variable Angle Locking Screws, self-tapping	Φ2.7×14	2
31477116	Variable Angle Locking Screws, self-tapping	Φ2.7×16	2
31477118	Variable Angle Locking Screws, self-tapping	Φ2.7×18	2
31477120	Variable Angle Locking Screws, self-tapping	Φ2.7×20	2
31477122	Variable Angle Locking Screws, self-tapping	Φ2.7×22	2
31477124	Variable Angle Locking Screws, self-tapping	Φ2.7×24	2
31477126	Variable Angle Locking Screws, self-tapping	Φ2.7×26	2
31477128	Variable Angle Locking Screws, self-tapping	Φ2.7×28	2
31477130	Variable Angle Locking Screws, self-tapping	Φ2.7×30	2
31477132	Variable Angle Locking Screws, self-tapping	Φ2.7×32	2
31477134	Variable Angle Locking Screws, self-tapping	Φ2.7×34	2
31477136	Variable Angle Locking Screws, self-tapping	Φ2.7×36	2
31477138	Variable Angle Locking Screws, self-tapping	Φ2.7×38	2
31477140	Variable Angle Locking Screws, self-tapping	Φ2.7×40	2
31478110	Variable Angle Locking Screws, self-tapping	Φ3.5×10	2
31478112	Variable Angle Locking Screws, self-tapping	Φ3.5×12	2
31478114	Variable Angle Locking Screws, self-tapping	Φ3.5×14	2
31478116	Variable Angle Locking Screws, self-tapping	Φ3.5×16	2
31478118	Variable Angle Locking Screws, self-tapping	Φ3.5×18	2
31478120	Variable Angle Locking Screws, self-tapping	Φ3.5×20	2
31478122	Variable Angle Locking Screws, self-tapping	Φ3.5×22	2
31478124	Variable Angle Locking Screws, self-tapping	Φ3.5×24	2
31478126	Variable Angle Locking Screws, self-tapping	Φ3.5×26	2
31478128	Variable Angle Locking Screws, self-tapping	Φ3.5×28	2
31478130	Variable Angle Locking Screws, self-tapping	Φ3.5×30	2
31478132	Variable Angle Locking Screws, self-tapping	Φ3.5×32	2
31478134	Variable Angle Locking Screws, self-tapping	Φ3.5×34	2
31478136	Variable Angle Locking Screws, self-tapping	Φ3.5×36	2
31478138	Variable Angle Locking Screws, self-tapping	Φ3.5×38	2
31478140	Variable Angle Locking Screws, self-tapping	Φ3.5×40	2
31478146	Variable Angle Locking Screws, self-tapping	Φ3.5×46	2
31478150	Variable Angle Locking Screws, self-tapping	Φ3.5×50	2
31478156	Variable Angle Locking Screws, self-tapping	Φ3.5×56	2

INSTRUMENTS

SCREW CASE

Code	Description	Size	QTY.
31478160	Variable Angle Locking Screws, self-tapping	Φ3.5×60	2
31478166	Variable Angle Locking Screws, self-tapping	Φ3.5×66	2
31478170	Variable Angle Locking Screws, self-tapping	Φ3.5×70	2
31478176	Variable Angle Locking Screws, self-tapping	Φ3.5×76	2
31478180	Variable Angle Locking Screws, self-tapping	Φ3.5×80	2
31475232	Cannulated Locking Screw, partial threaded	Φ4.0×32	2
31475240	Cannulated Locking Screw, partial threaded	Φ4.0×40	2
31475248	Cannulated Locking Screw, partial threaded	Φ4.0×48	2
31475255	Cannulated Locking Screw, partial threaded	Φ4.0×55	2
31475258	Cannulated Locking Screw, partial threaded	Φ4.0×58	2
31475262	Cannulated Locking Screw, partial threaded	Φ4.0×62	2
31475265	Cannulated Locking Screw, partial threaded	Φ4.0×65	2
31475270	Cannulated Locking Screw, partial threaded	Φ4.0×70	2
31475272	Cannulated Locking Screw, partial threaded	Φ4.0×72	2
31475278	Cannulated Locking Screw, partial threaded	Φ4.0×78	2
31475280	Cannulated Locking Screw, partial threaded	Φ4.0×80	2
31475285	Cannulated Locking Screw, partial threaded	Φ4.0×85	2
31475288	Cannulated Locking Screw, partial threaded	Φ4.0×88	2

INSTRUMENTS

SCREW CASE

291050 VariLoc 5.0mm Screw Case



Code	Description	Size	QTY.
30386012	Cortex Screw, fully threaded,self-tapping	Φ4.5×12	3
30386014	Cortex Screw, fully threaded,self-tapping	Φ4.5×14	3
30386016	Cortex Screw, fully threaded,self-tapping	Φ4.5×16	3
30386018	Cortex Screw, fully threaded,self-tapping	Φ4.5×18	3
30386020	Cortex Screw, fully threaded,self-tapping	Φ4.5×20	3
30386022	Cortex Screw, fully threaded,self-tapping	Φ4.5×22	3
30386024	Cortex Screw, fully threaded,self-tapping	Φ4.5×24	3
30386026	Cortex Screw, fully threaded,self-tapping	Φ4.5×26	3
30386028	Cortex Screw, fully threaded,self-tapping	Φ4.5×28	3
30386030	Cortex Screw, fully threaded,self-tapping	Φ4.5×30	3
30386032	Cortex Screw, fully threaded,self-tapping	Φ4.5×32	3
30386034	Cortex Screw, fully threaded,self-tapping	Φ4.5×34	3
30386036	Cortex Screw, fully threaded,self-tapping	Φ4.5×36	3
30386038	Cortex Screw, fully threaded,self-tapping	Φ4.5×38	3
30386040	Cortex Screw, fully threaded,self-tapping	Φ4.5×40	3
30386042	Cortex Screw, fully threaded,self-tapping	Φ4.5×42	3
30386044	Cortex Screw, fully threaded,self-tapping	Φ4.5×44	3
30386046	Cortex Screw, fully threaded,self-tapping	Φ4.5×46	3
30386048	Cortex Screw, fully threaded,self-tapping	Φ4.5×48	3
30386050	Cortex Screw, fully threaded,self-tapping	Φ4.5×50	3
30386052	Cortex Screw, fully threaded,self-tapping	Φ4.5×52	3
30386054	Cortex Screw, fully threaded,self-tapping	Φ4.5×54	3
30386056	Cortex Screw, fully threaded,self-tapping	Φ4.5×56	3
30386058	Cortex Screw, fully threaded,self-tapping	Φ4.5×58	3
30386060	Cortex Screw, fully threaded,self-tapping	Φ4.5×60	3
30386062	Cortex Screw, fully threaded,self-tapping	Φ4.5×62	3
30386065	Cortex Screw, fully threaded,self-tapping	Φ4.5×65	3
30386070	Cortex Screw, fully threaded,self-tapping	Φ4.5×70	3
30386075	Cortex Screw, fully threaded,self-tapping	Φ4.5×75	3
30386080	Cortex Screw, fully threaded,self-tapping	Φ4.5×80	3
34255014	Locking Screw, self-tapping	Φ5.0×14	6
34255016	Locking Screw, self-tapping	Φ5.0×16	6
34255018	Locking Screw, self-tapping	Φ5.0×18	6
34255020	Locking Screw, self-tapping	Φ5.0×20	6
34255022	Locking Screw, self-tapping	Φ5.0×22	6
34255024	Locking Screw, self-tapping	Φ5.0×24	6
34255026	Locking Screw, self-tapping	Φ5.0×26	6
34255028	Locking Screw, self-tapping	Φ5.0×28	6
34255030	Locking Screw, self-tapping	Φ5.0×30	6
34255032	Locking Screw, self-tapping	Φ5.0×32	6
34255034	Locking Screw, self-tapping	Φ5.0×34	6

INSTRUMENTS

SCREW CASE

Code	Description	Size	QTY.
34255036	Locking Screw, self-tapping	Ø5.0×36	6
34255038	Locking Screw, self-tapping	Ø5.0×38	6
34255040	Locking Screw, self-tapping	Ø5.0×40	6
34255042	Locking Screw, self-tapping	Ø5.0×42	6
34255044	Locking Screw, self-tapping	Ø5.0×44	6
34255046	Locking Screw, self-tapping	Ø5.0×46	6
34255048	Locking Screw, self-tapping	Ø5.0×48	6
34255050	Locking Screw, self-tapping	Ø5.0×50	6
34255055	Locking Screw, self-tapping	Ø5.0×55	6
34255060	Locking Screw, self-tapping	Ø5.0×60	6
34255065	Locking Screw, self-tapping	Ø5.0×65	6
34255070	Locking Screw, self-tapping	Ø5.0×70	6
34255075	Locking Screw, self-tapping	Ø5.0×75	6
34255080	Locking Screw, self-tapping	Ø5.0×80	6
34255085	Locking Screw, self-tapping	Ø5.0×85	6
34255090	Locking Screw, self-tapping	Ø5.0×90	6
30405120	Cancellous Screws, fully threaded	6.5x20mm	2
30405125	Cancellous Screws, fully threaded	6.5×25mm	2
30405130	Cancellous Screws, fully threaded	6.5×30mm	2
30405135	Cancellous Screws, fully threaded	6.5×35mm	2
30405140	Cancellous Screws, fully threaded	6.5×40mm	2
30405145	Cancellous Screws, fully threaded	6.5×45mm	2
30405150	Cancellous Screws, fully threaded	6.5×50mm	2
30405155	Cancellous Screws, fully threaded	6.5×55mm	2
30405160	Cancellous Screws, fully threaded	6.5×60mm	2
30405165	Cancellous Screws, fully threaded	6.5×65mm	2
30405170	Cancellous Screws, fully threaded	6.5×70mm	2
30405175	Cancellous Screws, fully threaded	6.5×75mm	2
30405180	Cancellous Screws, fully threaded	6.5×80mm	2
30405185	Cancellous Screws, fully threaded	6.5×85mm	2
30405190	Cancellous Screws, fully threaded	6.5×90mm	2
30405240	Cancellous Screws, partially threaded, 16mm Thread	6.5×40mm	2
30405245	Cancellous Screws, partially threaded, 16mm Thread	6.5×45mm	2
30405250	Cancellous Screws, partially threaded, 16mm Thread	6.5×50mm	2
30405255	Cancellous Screws, partially threaded, 16mm Thread	6.5×55mm	2
30405260	Cancellous Screws, partially threaded, 16mm Thread	6.5×60mm	2
30405265	Cancellous Screws, partially threaded, 16mm Thread	6.5×65mm	2
30405270	Cancellous Screws, partially threaded, 16mm Thread	6.5×70mm	2
30405275	Cancellous Screws, partially threaded, 16mm Thread	6.5×75mm	2
30405280	Cancellous Screws, partially threaded, 16mm Thread	6.5×80mm	2
30405285	Cancellous Screws, partially threaded, 16mm Thread	6.5×85mm	2

INSTRUMENTS

SCREW CASE

Code	Description	Size	QTY.
30405290	Cancellous Screws, partially threaded, 16mm Thread	6.5×90mm	2
31479112	Variable Angle Locking Screws, self-tapping	Ø5.0×12	2
31479115	Variable Angle Locking Screws, self-tapping	Ø5.0×15	2
31479118	Variable Angle Locking Screws, self-tapping	Ø5.0×18	2
31479120	Variable Angle Locking Screws, self-tapping	Ø5.0×20	2
31479122	Variable Angle Locking Screws, self-tapping	Ø5.0×22	2
31479125	Variable Angle Locking Screws, self-tapping	Ø5.0×25	2
31479128	Variable Angle Locking Screws, self-tapping	Ø5.0×28	2
31479130	Variable Angle Locking Screws, self-tapping	Ø5.0×30	2
31479132	Variable Angle Locking Screws, self-tapping	Ø5.0×32	2
31479135	Variable Angle Locking Screws, self-tapping	Ø5.0×35	2
31479138	Variable Angle Locking Screws, self-tapping	Ø5.0×38	2
31479140	Variable Angle Locking Screws, self-tapping	Ø5.0×40	2
31479142	Variable Angle Locking Screws, self-tapping	Ø5.0×42	2
31479145	Variable Angle Locking Screws, self-tapping	Ø5.0×45	2
31479148	Variable Angle Locking Screws, self-tapping	Ø5.0×48	2
31479150	Variable Angle Locking Screws, self-tapping	Ø5.0×50	2
31479155	Variable Angle Locking Screws, self-tapping	Ø5.0×55	2
31479160	Variable Angle Locking Screws, self-tapping	Ø5.0×60	2
31479165	Variable Angle Locking Screws, self-tapping	Ø5.0×65	2
31479170	Variable Angle Locking Screws, self-tapping	Ø5.0×70	2
31479175	Variable Angle Locking Screws, self-tapping	Ø5.0×75	2
31479180	Variable Angle Locking Screws, self-tapping	Ø5.0×80	2
31479185	Variable Angle Locking Screws, self-tapping	Ø5.0×85	2
31479190	Variable Angle Locking Screws, self-tapping	Ø5.0×90	2
31479195	Variable Angle Locking Screws, self-tapping	Ø5.0×95	2
31479196	Variable Angle Locking Screws, self-tapping	Ø5.0×100	2
31476250	Cannulated Locking Screw, partial threaded	Ø5.5×50	2
31476255	Cannulated Locking Screw, partial threaded	Ø5.5×55	2
31476262	Cannulated Locking Screw, partial threaded	Ø5.5×62	2
31476265	Cannulated Locking Screw, partial threaded	Ø5.5×65	2
31476270	Cannulated Locking Screw, partial threaded	Ø5.5×70	2
31476272	Cannulated Locking Screw, partial threaded	Ø5.5×72	2
31476278	Cannulated Locking Screw, partial threaded	Ø5.5×78	2
31476280	Cannulated Locking Screw, partial threaded	Ø5.5×80	2
31476285	Cannulated Locking Screw, partial threaded	Ø5.5×85	2
31476292	Cannulated Locking Screw, partial threaded	Ø5.5×92	2

INSTRUMENTS

SCREW CASE

293050 5.0/7.0mm Screw Case



Code	Description	Size	QTY.
34255030	Locking Screw, self-tapping	Φ5.0×30	2
34255032	Locking Screw, self-tapping	Φ5.0×32	2
34255034	Locking Screw, self-tapping	Φ5.0×34	2
34255035	Locking Screw, self-tapping	Φ5.0×35	2
34255036	Locking Screw, self-tapping	Φ5.0×36	2
34255038	Locking Screw, self-tapping	Φ5.0×38	2
34255040	Locking Screw, self-tapping	Φ5.0×40	2
34255042	Locking Screw, self-tapping	Φ5.0×42	2
34255044	Locking Screw, self-tapping	Φ5.0×44	2
34255045	Locking Screw, self-tapping	Φ5.0×45	2
34255046	Locking Screw, self-tapping	Φ5.0×46	2
34255048	Locking Screw, self-tapping	Φ5.0×48	2
34255050	Locking Screw, self-tapping	Φ5.0×50	2
31474070	Cannulated Locking Screw, partial threaded	Φ7.0×70	3
31474075	Cannulated Locking Screw, partial threaded	Φ7.0×75	3
31474078	Cannulated Locking Screw, partial threaded	Φ7.0×78	3
31474082	Cannulated Locking Screw, partial threaded	Φ7.0×82	3
31474085	Cannulated Locking Screw, partial threaded	Φ7.0×85	3
31474090	Cannulated Locking Screw, partial threaded	Φ7.0×90	3
31474095	Cannulated Locking Screw, partial threaded	Φ7.0×95	3
31474100	Cannulated Locking Screw, partial threaded	Φ7.0×100	3
31474105	Cannulated Locking Screw, partial threaded	Φ7.0×105	3
31474110	Cannulated Locking Screw, partial threaded	Φ7.0×110	3
31474115	Cannulated Locking Screw, partial threaded	Φ7.0×115	3
31474120	Cannulated Locking Screw, partial threaded	Φ7.0×120	3