

LCP Clavicular Reconstruction Plate


Product No.	Hole	Direction	Length
22563106	6	Right	102
22563107	7	Right	112
22563108	8	Right	121
22563206	6	Left	102
22563207	7	Left	112
22563208	8	Left	121

LCP Clavicle Hook

Product No.	Hole	Direction	Length
22586102	3+2	Right	45
22586103	3+3	Right	55
22586105	3+5	Right	75
22586107	3+7	Right	95
22586202	3+2	Left	45
22586203	3+3	Left	55
22586205	3+5	Left	75
22586207	3+7	Left	95

LCP Anterolateral Clavicular Plate

Product No.	Hole	Direction	Length
22578110	10	Right	68
22578111	11	Right	80
22578112	12	Right	92
22578113	13	Right	104
22578114	14	Right	116
22578115	15	Right	128
22578210	10	Left	68
22578211	11	Left	80
22578212	12	Left	92
22578213	13	Left	104
22578214	14	Left	116
22578215	15	Left	128



Monoloc
LCP Clavicle Plates

CHANGZHOU KANGHUI MEDICAL INNOVATION CO., LTD

Add: No.11 North Changjiang Road, Xinbei Zone, Changzhou, Jiangsu 213022, P.R.China
Tel: +86-519-85139851 Fax: +86-519-85128628
www.kanghui.com
Shanghai Office
Add: Building 18, No. 1000, Jinhai Road, PuDong District, Shanghai, China 201206
Tel: +86-21-50319916 Fax: +86-21-50312913

Monoloc

LCP Clavicle Plates

Indications



LCP Superior Anterior Clavicle Plate without lateral extension

- Fractures of the clavicle shaft
- Malunions of the clavicle shaft
- Non-unions of the clavicle shaft



LCP Superior Anterior Clavicle Plate with lateral extension

- Fractures of the lateral clavicle
- Malunions of the lateral clavicle
- Non-unions of the lateral clavicle
- Shaft fractures



LCP Clavicle Hook Plate

- Acromioclavicular joint dislocation
- Fractures of the lateral clavicle

Key Features & Benefits

- Anatomically pre-shaped, twisted die sign
- Combines strength with ductility
- Left and right plates
- Notches in the plate allow additional plate contouring
- Limited-contact shaft profile
- Tapered tip for sub muscular plate insertion
- Plates with lateral extension feature six distal locking holes that accept 2.7mm locking and 2.4mm cortex screws
- The holes in the plate shaft accept 3.5mm locking, 3.5mm cortex, and 4.0mm cancellous bone screws

