


# VariLoc

Solutions to Clavicle Fractures




THERAPY OPTIMIZED  
**EXPECT EXTRAORDINARY**

## Clavicular Locking Compression Plate System


### Indication

Variloc Clavicle Locking Plate is indicated for fixation of fractures, malunions, nonunions and osteotomies of the clavicle, as well as acromioclavicular joint injuries.



These images do not represent product efficacy or final result and should be used as a pictorial representation only.

### Unique Design: Variable Angle Locking Screw




10°-2.4/2.7/3.5/5.0mm VA screw

**This enables the surgeon to:**

- Adapt screw angulations to patient anatomy
- Capture specific fracture fragments
- Adjust screw trajectory after plate bending
- Position screws precisely to avoid joint penetration

### Features & Benefits

- Recon design allows for plate contouring to fit patient anatomy
- Percutaneous insertion prevents soft issue irritation
- Undercuts reduce impairment of blood supply
- Anatomically pre-contoured shape facilitates optimal implant placement and surgery
- Rounded shaft profile minimizes the risk of conflicts between the plate and surrounding soft tissue
- Intra-operative choice of clavicle hook size optimizes size matching for each individual patient



3 different hook depths

18mm  
15mm  
12mm

