


VariLoc

Solutions to Femoral Neck Fractures



THERAPY OPTIMIZED
EXPECT EXTRAORDINARY

Proximal Femoral Locking Plate with 7.0mm Cannulated Locking Screws




Locking Plate

- Designed to secure optimal fit to anatomical shape of the femoral lateral wall
- Minimum invasive and less blood loss
- Left and right option to secure fit to the femoral anteversion
- Positioned on the flat lateral wall beneath the greater trochanter for easy positioning
- The 130° CCD angle fits to the anatomical angle to minimize the risk of screw breakage

7.0mm Cannulated Locking Screws

- Prevent screws pull out and femoral head dispositioning
- Reduce the risk of screwdriver breakage
- Enhance effective compression



5.0mm Locking Screws

- Secure strong shaft fixation
- Reduce the risk of screw driver breakage
- Enhance the screw pull out strength